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
Quick Reference Card

for dental practice

Orientation guidelines for the application of BioHPP

- Preparing & modelling
- Sterilising, conditioning & fitting
- Cleaning at the practice
- Cleaning in the laboratory
- Cleaning at the patient's home


Key

visio.lign®
The Aesthetic and Functional system

BioHPP

Tooth structure/
implant

Preparing & modelling

Sterilising, conditioning
& fitting

Cleaning at the practice

Cleaning in the laboratory

Cleaning at the patient's
home

Very well suited	Suitable	Possible	Unsuitable
✓✓	✓	○	✗

In the
dental practice


In the
dental laboratory

At the
patient's home

Preparation

	Chamfer preparation	Step preparation	Tangential preparation
BioHPP	✓	✓	✗
HIPC	✓	✓	✗


Preparation


Intraoral or extraoral preparation? _____

Flap: Extra-oral preparation of the abutment

Flapless: Intra-oral preparation of the abutment

_____ Intraoral or
extraoral
preparation? _____

Instruments

for the preparation of BioHPP elegance abutments


Intra-oral:

Composite
Q-finishing bur
from Komet:
H375RQ.314.016
H379Q.314.018


Use with spray
cooling (at least
50 ml/min)
20,000 rpm

Extra-oral:

Milling set
BioHPP Perfor-
mance
Polymers® 7-pc

REF 33000831

Use of a laboratory analogue.
Horizontal grooves serve as
mechanical retention during
veneering


Instruments

for the preparation of BioHPP
elegance abutments

Temporary restoration

BioHPP elegance Abutment

In the temporary restoration only the height of the abutment is reduced.
The sulcus course will only be customised chairside in the case of permanent restoration.



Temporary restoration

BioHPP elegance Abutment

Modelling

for temporary and permanent restoration



BioHPP can be scanned intraorally without powder.

Modelling

Sterilising & disinfecting

Tooth structure/
implant

	Sterilising (autoclaving, max. 160°C)	Disinfecting (e.g. chlorhexidine, chlorhexamed, alcohol)
BioHPP elegance abutment only	✓	✓
BioHPP as a veneered crown abutment	✗	✓
BioHPP bonded to titanium adhesive bases (not recommended by bredent)	✗	✓


Sterilising &
disinfecting

Conditioning & fitting

temporary

Crown-based approach

BioHPP crown on various implant materials or tooth

Conditioning & fitting

temporary

Crown-based approach ·
BioHPP crown on various
framework materials

Step	1	2	3				
	BioHPP Blast crown with 110 µm	visio.link on BioHPP	Bonding	Metal/ alloys	Zirconium dioxide	BioHPP	Dentine/ enamel
		X	Eugenol-free zinc oxide cement (Tempbond (Kerr))	✓	✓	✓	○*
		X	A-Silicone-based fixing cement (Temposil 2 (Coltene Whaledent))	✓✓	✓✓	✓✓	✓✓


* only use with a preparation angle of up to 5°

Conditioning & fitting

permanent

Crown-based approach

BioHPP crown on various teeth/implant/metal alloys

Conditioning & fitting

permanent

Crown-based approach ·
BioHPP crown on various
framework materials

Step

1

2

3

BioHPP Blast crown with 110 µm	visio.link on BioHPP	Bonding	Metal/ alloys	Zirconium dioxide	BioHPP	Dentine/ enamel
	✓✓	Adhesive technique (VarioLink; NX-3, Panavia 2.0 & V5)	✓✓	✓✓	✓✓	✓✓
	●	Self-adhesive (Rely X unicem (3M Espe))	✓	✓	✓	✗
✗		Glass ionomer cement (Ketac Cem (3M Espe))	* ●	● *	●	● *
✗		Zinc phosphate cement (Harvard Cement)	●	●	●	● *

* only use with a preparation angle of up to 5°

Conditioning & fitting

temporary

Framework-based approach

BioHPP framework together with different crown materials

Step	1	2	3					
	BioHPP Blast frame- work/abut- ment with 110 µm	visio.link on BioHPP	Bonding	Metal/ alloys	Zirco- nium dioxide	BioHPP	Lithium disilicate (e.g. emax)	HIPC
		X	Eugenol-free zinc oxide cement (Tempbond (Kerr))	✓*	✓*		X	✓
		X	A-Silicone-based fixing cement (Temposil 2 (Coltene Whaledent))	✓✓	✓✓	✓✓	X	✓

* only use with a preparation angle of up to 5°

Conditioning & fitting

temporary

Framework-based approach
BioHPP framework together
with different crown
materials

Conditioning & fitting

permanent

Framework-based approach


BioHPP framework together with different crown materials

Step

1

2

3

Blast BioHPP framework/ abutment with 110 µm	visio.link on BioHPP	Bonding	Metal/ alloys	Zirco- nium dioxide (1)	BioHPP	Lithium disilicate (e.g. emax)	HPC
	✓✓	Adhesive technique (Panavia 2.0 & V5 (Kuraray))	✓✓	✓✓	✓✓	✓✓ (K)	✓✓
	○	Self-adhesive (Rely X unicem (3M Espe))	✓	✓	✓	✗	○
	✗	Glass ionomer cement (Ketac Cem (3M Espe))	○*	○*	○	✗	○
	✗	Zinc phosphate cement (Harvard Cement)	○	○	○	✗	○

(1) Flexural strength > 900 MPa (K) only recommended for crowns * Only use for a preparation angle up to 5°


Conditioning & fitting

permanent

Framework-based approach
BioHPP framework together
with different crown
materials

Professional dental cleaning

subgingival* – fixed

	Scalers/ curettes	Ultrasonic scalars		Air-powder-water jets		
		Metal tips	Plastic tips/ carbon (e.g. Periosoft tips from Acteon)	Nozzle subgingival (e.g. PERIO-FLOW® nozzles from EMS or Perio Easy subgingival nozzle from Acteon)	Sodium carbonate powder (e.g. AIR- FLOW® Powder CLASSIC from EMS)	Glycine powder, approx. 25 µm grain size (e.g. AIR-FLOW® Powder PERIO from EMS or Air-N-Go "Perio" from Acteon)
BioHPP abutments	✗	✗	✗	✓	✗	✓
Titanium abutments	✗	✗	✓	✓	✗	✓
HIPC Kronenrand	✗	✗	✓	✓	✗	✓

* Cleaning study, Ms. Braun, 2015


Professional
dental cleaning

subgingival – fixed

Professional dental cleaning

supragingival – fixed

	Scalers/ curettes	Ultrasonic scalers		Air-powder-water jets			Polishing	
	Periosoft scaler (Kerr)	Metal tips	Plastic tips (e.g. KaVo Sonicflex Paro straight)	Nozzle supragingival (e.g. Supranozzle from Acteon)	Sodium carbonate: (e.g. Airflow comfort from EMS)	Erythritol powder (e.g. Airflow Plus from EMS)	Super- Snap (Shofu)*	Prisma gloss (Dentsply)*
BioHPP	✓	✗	✓	✓	✗	✓	✓✓	✓
HIPC/ novo.lign/ neo.lign	✓	✗	✓	✓	✗	✓✓	✓✓	✓
crea.lign	✓	✗	✓	✓	✓	✓✓	✓✓	✓

Sources:


Cleaning study in cooperation with the LMU Munich under the management of Dr. Bogna Stawarczyk

** The Journal of Prosthetic Dentistry: Surface properties of polyetheretherketone after different lab-and chairside polishing protocols

Professional
dental cleaning




supragingival – fixed

Cleaning in the laboratory

temporary

Removable dental prosthesis · (e.g. secondary telescope prosthesis, bar-based prosthesis)

	Needle cleaning unit (e.g. Sympro and Symprofluid)	Ultrasonic cleaning baths	Fine sandblasting unit (grain size: 50 µm)
BioHPP	✓✓	✓✓	 (Polishing required)
HIPC/ novo.lign/ neo.lign	✓✓	✓✓	 (Polishing required)
crea.lign	✓✓	✓✓	 (Polishing required)

Cleaning in
the laboratory

temporary


Removable dental
prosthesis · (e.g. secondary
telescope prosthesis,
bar-based prosthesis)

Cleaning at the patient's home

Fixed and removable restoration

Removable dental prosthesis

(e.g. secondary telescope prosthesis, bar-based prosthesis)

	Vibration cleaning baths	Soft toothbrushes	Medium toothbrushes	Firm toothbrushes	Electric toothbrushes	Ultrasonic toothbrushes	Cleaning rods
BioHPP	✓✓	✓✓	✓	✗	✓	✗	✓
HIPC/ novo.lign/ neo.lign	✓✓	✓✓	✓	✗	✓	✗	✓
crea.lign	✓✓	✓✓	✓	✗	✓	✗	✓

Cleaning

fixed and removable restoration

Removable dental prosthesis · (e.g. secondary telescope prosthesis, bar-based prosthesis)

Chipping Repair Kit

Chips can be restored easily and quickly

Chipping frequent in the case of ceramic veneers

Following repair with crea.lign® composite

Chipping Repair Kit

Primer works together with all prosthetic materials


crea.lign® incisal materials and dentine (in A2 and A3 for covering in 80% of cases)

100 treatments in one Chipping Repair Kit

No hydrofluoric acid necessary

The physiological solution

A restoration with BioHPP significantly reduces the chewing force peaks in both the vertical and the lateral direction of chewing movement. This cushioning property is both comfortable and natural for the patient and prolongs the life of the restoration.

Natural tooth.


A restoration with BioHPP alleviates the loading caused by natural and prosthetic-related forces on the implant.


Rigid materials prevent natural movement and transfer the force to the implants and the jaw bones.


Other literature that may be of interest to you ———


REF 000 547 GB

REF 000 534 GB

REF 000 577 GB

REF 000 651 GB

Other
literature that
may be
of interest
to you

Quick Reference Card

for dental practice

bredent
group

bredent GmbH & Co. KG
Weissenhorner Str. 2 | 89250 Senden | Germany
T: (+49) 0 73 09 / 8 72-4 41
F: (+49) 0 73 09 / 8 72-4 44
www.bredent.com · @: info@bredent.com

Mistake and subject
to change reserved
20170309