

PRINTER SETTINGS MASTER MODEL & CLASSIC MODEL / DRUCKEREINSTELLUNGEN MASTER MODEL & CLASSIC MODEL

This information serves as a guideline. Please make adjustments if necessary.
Diese Angaben dienen als Richtwert. Bei Bedarf bitte Anpassungen vornehmen.



Asiga UV Max 385	50 Micron	100 Micron	Curing	Settings
Layer Exposure (s)	0,91	1,13	30 min	Denta Model (Modified)
Base Exposure (s)	8,9	8,9		
Light Intensity mw/cm ²	11,6	11,6		
UV Dose mws/cm ²	10,60	13,1		



Phrozen Sonic Mini 4K	50 Micron	100 Micron	Curing	Settings
Layer Exposure (s)	5,5	5,5	30 min	Manual
Base Exposure (s)	40	40		

PRINTER SETTINGS MASTER MODEL & CLASSIC MODEL / DRUCKEREINSTELLUNGEN MASTER MODEL & CLASSIC MODEL

This information serves as a guideline. Please make adjustments if necessary.
Diese Angaben dienen als Richtwert. Bei Bedarf bitte Anpassungen vornehmen.



Anycubic Photon Mono 4K	50 Micron	100 Micron	Curing	Settings
Layer Exposure (s)	3,2	-	30 min	Manual
Base Exposure (s)	25	-		



Ackuretta SOL	70 Micron	100 Micron	Curing	Settings
Layer Exposure (s)	3	3,3	10 min	Curo Model / Nextdent Model 2.0
Base Exposure (s)	19,8	19,8		

PRINTER SETTINGS DENTAL CAST DRUCKEREINSTELLUNGEN DENTAL CAST

This information serves as a guideline. Please make adjustments if necessary.
Diese Angaben dienen als Richtwert. Bei Bedarf bitte Anpassungen vornehmen.



Asiga UV Max 385	50 Micron	100 Micron	Curing	Settings
Layer Exposure (s)	2,1	3,2	30 min	Nextdent Cast V2.2
Base Exposure (s)	8,8	8,8		
Light Intensity mw/cm ²	5	5		
UV Dose mws/cm ²	10,5	10,5		



Phrozen Sonic Mini 4K	50 Micron	100 Micron	Curing	Settings
Layer Exposure (s)	6	9	30 min	Manual
Base Exposure (s)	40	40		

PRINTER SETTINGS DENTAL CAST

DRUCKEREINSTELLUNGEN DENTAL CAST

This information serves as a guideline. Please make adjustments if necessary.
Diese Angaben dienen als Richtwert. Bei Bedarf bitte Anpassungen vornehmen.



Anycubic Photon Mono 4K	50 Micron	100 Micron	Curing	Settings
Layer Exposure (s)	3	-	30 min	Manual
Base Exposure (s)	30	-		



Ackuretta SOL	70 Micron	100 Micron	Curing	Settings
Layer Exposure (s)	2,4	-	10 min	Optiprint Cast / Whipmix-Vericast OS
Base Exposure (s)	11,2	-		