





- 1 The neck of the tooth is coated with the dentine mass A3.5 or beige modifier – one degree darker than the subsequent tooth colour.
- 2 The body of the tooth and mamelons are coated with the dentine mass A3.
- 3 Almost the entire edge is built up with Enamel E2.
- 4 The form of the cutting edge is completed with Incisal opal.

- 1 The neck of the tooth is coated with the dentine mass A3.5 or beige modifier – one degree darker than the subsequent tooth colour.
- 2 The body of the tooth and mamelons are coated with the dentine mass A3.
- 3 Add thin enamel cracks using visio.lign color brown.
- 4 Place alternating layers of Incisal blue and Incisal rose over the mamelons.
- 5 Apply Incisal blue in a mesial and distal direction from the cutting edge.
- 6 The entire edge is built up with Enamel E2.



Standard layering



Customised layering

Assignment of the crea.lign opaquer to dentine colours

crea.lign opaquer	1	2	3	4	5	6	7	8	9	Gum	Z
system shades	A1 B2	A2	A3	B1/C1/ BL3	C2/C3/ D2/D4	B3 B4	A3,5	A4 C4	D3	Gingiva shades	Zirconium liner

Assignment of the edges to dentine colours

Enamel / A-D Combination	BL3	A1	A2	A3	A3,5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
E1	X	X					X				X				X		
E2			X	X				X									
E3					X				X	X		X	X			X	X
E4						X								X			

Quality check on the veneer

- Is the layer of lubricant still present?
- Remove the layer of lubricant using crea.lign surface cleaner and a toothbrush.
- After steaming off the polishing residues, the veneer should have a high gloss.

Tips and tricks

Little space for the veneer:

- Combo.lign in the matching colour is applied as the first layer. Cover the combo.lign completely with the following layers; it must not lie on the surface, to prevent discolouration and plaque accumulation.

Sensitive patients (allergy patients):

- In order to ensure complete polymerisation of the veneer, the specified polymerisation times should be respected. Either place the completed work in the pressure pot for 1 hour at 55° C (no pressurisation required), or pour boiling water into a container and place the work in it for 30 minutes.



009837GB-20240228
Mistake and subject to change reserved

