



**visio.lign**

The Aesthetic and Functional System

3

# Presentation of the system

veneering composite



**crea.lign**

veneers



**novo.lign**

denture teeth



**neo.lign**

composite blanks



**visio.CAM**

stains



**visio.paint**

bonding system



**bond.lign**

perform with  
**visio.lign**





## Chapter overview

1

Preparatory work  
prepare natural beauty



silicones for matrices

08

2

crea.lign  
create natural beauty



veneering composite  
modelling liquid  
Sets

12

18

21

3

visio.paint  
colour natural beauty



visio.lign color  
visio.lign shield

28

28

4

novo.lign  
veneer natural beauty



Veneers  
Sets

34

43

5

neo.lign  
choose natural beauty



Teeth  
Sets

48

57

## 6

### bond.lign secure natural beauty



Primer and Bonder	60
Opaker	62
luting composite	64
DTK adhesive	64
Sets	66

## 7

### visio.CAM design natural beauty



Composite blanks	70
------------------	----

## 8

### Tools support natural beauty



bre.Lux	74
visio.lign Toolkit	78
crea.lign gel	79
crea.lign paste	80
novo.lign	80
bond.lign	81
Preparatory work	81
visio.lign tooth cabinet / box	82
Material processing matrix	83
Media	84
Index	86

1

# Preparatory work

prepare natural beauty

---



PROPERTIES SUCH AS PRECISE  
APPLICATION, DETAILED  
IMPRESSIONS AND LIGHT  
TRANSMISSION SPEAK IN  
FAVOUR OF THE FINE-GRAINED  
VISIO.SIL MATRIX SILICONES.

SILICONES FOR MATRICES

# Preparatory work

## silicones for matrices



### visio.sil fix

High-precision silicone for matrices  
Hardness of approx. 55 Shore A

visio.sil fix is an addition-curing, high-precision silicone for matrices. With visio.sil fix, the approximal spaces and gingival contouring of the wax-up are perfectly



Yellow

REF 54001300

50 ml


moulded to the finest detail. Once the restoration has been manufactured, it just has to be polished to a high gloss, saving effort and time.



### visio.sil

Transparent silicone for matrices  
Hardness of approx. 60 Shore A

visio.sil is a transparent, addition-curing silicone for matrices, used for the fabrication of translucent matrices and processing of light-curing materials.



transparent

REF 54001200

50 ml

#### Benefits of visio.sil and visio.sil ILT transparent silicones for matrices

visio.sil is a transparent, addition-curing silicone for matrices, used for the fabrication of translucent matrices and processing of light-curing materials.





## visio.sil ILT

Transparent silicone for matrices for Inverse Layering Technique

Hardness of approx. 75 Shore A

visio.sil ILT has been specially developed for the Inverse Layering Technique. This transparent, addition-curing silicone for matrices is characterised by its extreme hardness of approx. 75 Shore A. It is flowable so it fills the proximal space and

ensures perfect reproduction of the surface texture; however, it is also firm and stable and can be applied very precisely. visio.sil ILT can be applied in the patient's mouth (class 1 medical device) in the hardened state.



transparent

REF 54001400

50 ml



### Benefit of transparent

#### visio.sil ILT silicone for matrices

visio.sil ILT silicone for matrices has been approved for oral application and is perfectly suited for positioning and luting indirect veneers.

2

---



# crea.lign

create natural beauty



THE LIGHT-CURING CERAMIC  
COMPOSITE IMPRESSES WITH  
NATURAL OPALESCENCE  
AND THE BEST LIGHT  
TRANSMISSION -  
FOR NATURAL BEAUTY

VENEERING COMPOSITE  
STAINS  
MODELLING LIQUID

# crea.lign | veneering composite

create natural beauty - with the light-cured ceramic composite for natural beauty

crea.lign is a light-curing ceramic composite for the permanent veneering of metal, ceramic and polymer frameworks. The exceptional material properties adapt elasticity and hardness of crea.lign to the various substructure materials.

crea.lign contains only nano particles and no ground glass filler. It consists of 50% opalescent ceramic filler and a high-strength oligomer matrix.

The ceramic fillers used and the avoidance of hard glass fillers increase abrasion resistance and prevent plaque build-up, discoloration and embrittlement of the material.

Agglomeration effects (clumping) in the material are avoided and a homogeneous and dense surface is achieved with the particle size of approx. 40 µm and optimal polishing properties. Thanks to the dense surface, water absorption is low, which in turn ensures mechanical stability.



## crea.lign | Coding

### crea.lign gel:

	REF CLFN_____
dentine, e.g. A3:	REF CLFN <b>DA30</b>
opaque dentin, e.g. A (A3/A35):	REF CLFN <b>OPDA335</b>
transpa clear, e.g. T1:	REF CLFN <b>00T1</b>
Incisal, e.g. I1:	REF CLFN <b>00I1</b>
Incisal Enamel, e.g. E1:	REF CLFN <b>00E1</b>
High Translucent, e.g. HT1 (A1/A2):	REF CLFN <b>HT01</b>
GUM, e.g. G1:	REF CLFN <b>00G1</b>
Modifier, e.g. M1:	REF CLFN <b>00M1</b>
Stains, e.g. S1:	REF CLFN <b>00S1</b>

### crea.lign paste:

	REF CLPND_____
dentine, e.g. A3:	REF CLPND <b>A30</b>
transpa opal, e.g. T1:	REF CLPN <b>00T1</b>
Incisal, e.g. E1:	REF CLPN <b>00E1</b>
GUM, e.g. G50:	REF CLPN <b>0G50</b>
–	–
–	–



## Natural beauty

- Natural opalescence
- Optimal depth effect thanks to coordinated refractive indices
- Translucency (greater degree of freedom than in ceramics)

crea.lign is available in both dosage forms gel and paste in classic A-D colours and other GUM and incisal materials and can be seamlessly combined with each other.

## crea.lign | Benefits

**Natural** - Balance between opalescence, fluorescence and translucency creates a natural, vibrant colour effect

**Physiology** - crea.lign is significantly more impact-resistant than ceramics thanks to its excellent mechanical properties. This increases the resilience of the veneer.

**Safety** - crea.lign's permanent chemical adhesion to all framework materials as well as plaque and discoloration resistance ensure the long-term success of the veneer.

**crea.lign gel** is the "more fluid" variant of the two dosage forms of crea.lign veneering composite. The specially developed single-hand grip and application cannula make precise dosing with just one hand easier for the user.

**crea.lign paste** is the pasty variant of crea.lign's two dosage forms. The consistency of crea.lign pastes is soft and smooth. Contouring is very easy. No undesired adhesion. The composite shaper kit makes modelling even easier. The silicone modelling instruments have been specially developed for the crea.lign paste and serve as an attachment and modelling aid. This makes it easy to achieve the finest transitions in modelling.

- All colours of the crea.lign gel or crea.lign paste can be found on page 16 and 17
- crea.lign opaque, see page 63
- shade combination tables, see page 18



crea.lign gel



**crea.lign** paste

# crea.lign | gel

BL3	A1	A2	A3	A3.5	A4	B1	B2	B3
<u>BL3</u>	<u>A10</u>	<u>A20</u>	<u>A30</u>	<u>A35</u>	<u>A40</u>	<u>B10</u>	<u>B20</u>	<u>B30</u>
	B4	C1	C2	C3	C4	D2	D3	D4
	<u>B40</u>	<u>C10</u>	<u>C20</u>	<u>C30</u>	<u>C40</u>	<u>D20</u>	<u>D30</u>	<u>D40</u>

A3/A3.5	B3	C3	D3
<u>A</u>	<u>B</u>	<u>C</u>	<u>D</u>
<u>A335</u>	<u>B3</u>	<u>C3</u>	<u>D3</u>
Bleach	Bleach (A1/B1)		
<u>B</u>	<u>L</u>	<u>A</u>	<u>B</u>
<u>1</u>			

E1	E2	E3	E4
<u>E1</u>	<u>E2</u>	<u>E3</u>	<u>E4</u>

opal	blue	rose	universal
<u>I1</u>	<u>I2</u>	<u>I3</u>	<u>I4</u>
<u>I5*</u>	*C.G.S. finish		
<u>I5</u>			

clear
<u>T1</u>
<u>T1</u>

A1/A2	A2/A3	A3/A3.5
<u>HT1</u>	<u>HT2</u>	<u>HT3</u>
<u>01</u>	<u>02</u>	<u>03</u>

light	pink	pink	red
<u>G1</u>	<u>G2</u>	<u>G3</u>	<u>G4</u>
<u>G1</u>	<u>G2</u>	<u>G3</u>	<u>G4</u>

REF

**CLFND** \_ \_



### gel Dentine

is available in the classic A-D dentine shades and BL3 bleached shade.

REF

**CLFNOPD** \_ \_ \_ \_



### Opaque Dentin Gel

crea.lign Opaque Dentine materials provide chroma and depth effect in composite veneers with low wall thicknesses. With the materials OD A to OD D and OD bleach the most common colors are made available.

REF

**CLFN00** \_ \_



### gel incisal

Alongside the enamel materials (E1/E2/E3/E4), crea.lign incisal (opal/blue/pink/universal/C.G.S.finish) can be used for specific accentuation of the incisal area (see also colour combination table of enamel and dentine materials on page 17)

REF

**CLFN00** \_ \_



### gel Transpa clear

crea.lign gel Transpa clear is the crystal-clear version of the crea.lign gel.

REF

**CLFNHT** \_ \_



### High Translucent Gel

crea.lign HT are chromatically attenuated Dentine materials with a very high degree of translucency. They allow a natural transition from the dentine to the enamel area and they are available in the versions HT1, HT2 and HT3.

REF

**CLFN00** \_ \_



### gel GUM

crea.lign gel GUM (light/pink/pink/red) are customisation materials for the gum area.





REF  
CLFN00 \_\_



gel **Modifier**

(beige/olive/caramel/purple/umbra) are intensively coloured materials, for example, to accentuate the neck region of a tooth or replicate decalcification stains.




REF  
CLFN00 \_\_



**Stains**

crea.lign stains (brown/orange) are light-curing intensive materials and are used for colour customisation of composite veneers. crea.lign stains must not be coated; they are used to emphasise, for example, proximal spaces, incisal edges and fissures.

# crea.lign | paste



REF  
CLPND \_\_



paste **dentine**

is available in the classic A-D dentine shades and BL3 bleached shade.



REF  
CLPN00 \_\_



paste **Incisal**

like crea.lign gel, is available in 4 shades (E1/E2/E3/E4) and is used to design the cutting area. The incisal materials are perfectly matched to the dentine materials (see table on this page below).



REF  
CLPN00 \_\_



paste **Transpa opal**

is the transparent variant of crea.lign paste. Thanks to optimal translucency and opalescence, it gives a natural, vibrant colour effect.



REF  
CLPNO \_\_



paste **GUM**

crea.lign paste GUM (light/opaque/bright/orange/brown) is used to coat the gingival elements for red-white customisation.

## crea.lign | shade combination tables

crea.lign gel/paste*	classic A-D dentine shades																
	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4	B13
Enamel																	
E1	■					■				■				■			■
E2		■	■				■										
E3				■				■	■		■	■			■	■	
E4					■								■				

\* All crea.lign gel materials can be seamlessly combined with all crea.lign paste materials.

## crea.lign | modelling liquid



During modelling instruments can be wetted with the crea.lign Modelling Liquid in order to avoid an adhesion of the veneering composite.

REF CLFMOD10

10 ml



**crea.lign** veneering system  
paste + gel,  
complete set



# crea.lign | sets

**REF CLIGNSETFK** crea.lign veneering system paste + gel, complete set

aluminium case,  
width 42.5 cm x height 43.5 cm x depth 20.5 cm



1 x

modelling liquid

10 ml

1 x



visio.link

10 ml

1 x



Disposable trays

1 x 12



Mixing pad

2 x



MKZ primer

5 ml

1 x



MKZ EM activator

4 ml

1 x



K primer

4 ml

1 x



Spatula

1 x



Assortment contains all available shades of the following components:

- crea.lign opaque
- crea.lign paste
  - dentine · Transpa opal
  - incisal · GUM
- crea.lign gel
  - incisal · Transpa clear
  - dentine · Modifier
  - GUM Opaque Dentin
  - High Translucent
- crea.lign Stain

➔ see page 16/17


Primer and activator

➔ see page 60/61

crea.lign syringe holder

stainless steel

2 x



Disposable brush

20 x




Brush holder

1 x




crea.lign paste incisal, 3 g

4 x




crea.lign paste dentine, 3 g

17 x




Rotary piston

27 x




crea.lign paste Transpa opal, 3 g

1 x




crea.lign paste GUM, 3 g

5 x



Single-hand grip crea.lign gel with UV LED

12 x



Pin cover

2 x 10




Application cannula

2 x 10




crea.lign gel incisal, 5 g

9 x




crea.lign Opaque Dentin, 5 g

5 x




crea.lign shade guide dentine shade (empty)

1 x




crea.lign Stain, 1,4 g

2 x




crea.lign gel Transpa clear, 5 g

1 x



crea.lign High Translucent, 5 g

3 x




crea.lign gel Modifier, 5 g

5 x




crea.lign opaque, 4 g

11 x




crea.lign gel dentine, 5 g

17 x




crea.lign gel GUM, 5 g

4 x



crea.lign shade guide Incisal, Stains & GUM (empty)


1 x



visio.lign information folder

with various printed materials, such as crea.lign freestyle manual and table display with illustrated layering instructions.

1 x



crea.lign Composite shaper kit

1 x angle chisel

1 x taper point

1 x flat chisel


1 x cup chisel

1 x cup round

**REF CLIGNSET12**

**crea.lign Assortment with syringe holder**

crea.lign syringe holder  
plastic



1 x

Single-hand grip crea.lign gel  
with UV LED



12 x

Pin cover  
2 x 10




Application cannula  
2 x 10




crea.lign Stain, 1,4 g

2 x



crea.lign gel, 5 g

12 x



crea.lign stains and gel  
shades of your choice  
➔ see page 16/17

modelling liquid, 10 ml

1 x




**REF CLPNSET9**

**crea.lign paste starter kit, 9 pastes can be selected**

crea.lign paste, 3 g

9 x



Rotary piston

9 x



crea.lign Composite shaper kit

1 x angle chisel

1 x laper point

1 x flat chisel

1 x cup chisel


1 x cup round

crea.lign paste  
shades of your choice  
➔ see page 17

**REF CLPNSET10**

**crea.lign paste starter kit, 10 x crea.lign paste**

syringe holder for crea.lign paste



2 x

crea.lign paste dentine, 3 g



- 1 x A2
- 1 x A3
- 1 x A3,5
- 1 x B3
- 1 x C3
- 1 x D3

crea.lign paste Enamel, 3 g



- 1 x E1
- 1 x E2
- 1 x E3

crea.lign paste Transpa opal,, 3 g



- 1 x T1

Mixing pad

1 x



Rotary piston




10 x

**REF CLPNSETFK**

**crea.lign paste complete kit**


crea.lign paste dentine, 3 g  
  
 17 x

Rotary piston  
  
 27 x


visio.link  
 10 ml  
  
 1 x


crea.lign paste shades of your choice  
 → see page 17


crea.lign paste GUM, 3 g  
  
 5 x


modelling liquid  
 10 ml  
  
 1 x

crea.lign paste incisal, 3 g  
  
 4 x




crea.lign opaque, 4 g  
  
 10 x



crea.lign surface cleaner  
 125 ml  
  
 1 x

crea.lign paste Transpa opal,, 3 g  
  
 1 x

crea.lign opaque Z,  
 zirconium liner 4 g  
  
 1 x

crea.lign Composite shaper kit

1 x angle chisel  
  
 1 x taper point  
  
 1 x flat chisel  


1 x cup chisel  
  
 1 x cup round  


**crea.lign freestyle intro Kit**

**REF CLFKSETA20**

**REF CLFKSETA30**


**REF CLFKSETA35**

A2  
 oder  
 A3  
 oder  
 A3,5

1 x crea.lign gel dentine, 5 g  


1 x crea.lign paste dentine, 3 g  


1 x crea.lign gel incisal, 5 g  


1 x crea.lign opaque, 4 g  


crea.lign gel and crea.lign opaque the set supplied corresponds to the assortment shade selected, A3 or A3.5)

**REF CLIGNSETG**


**Red-white aesthetics set**

visio.lign folding box



1 x

Layering instructions



1 x

Unique brush no. 4



1 x

modelling liquid  
10 ml



1 x

visio.link  
10 ml



1 x

Demo prosthesis



2 x

Pin cover

1 x 10



Application cannula

1 x 10



Mixing pad




1 x

Rotary piston



5 x

Single-hand grip crea.lign gel  
with UV LED



1 x


crea.lign gel incisal, 5 g



1 x I1 opal

1 x I2 blue

crea.lign gel Transpa clear, 5 g



1 x T1 transpa clear


crea.lign gel Modifier, 5 g



1 x M1 beige

1 x M4 purple

crea.lign gel GUM, 5 g




1 x G1 light

1 x G2 pink

1 x G3 pink


1 x G4 red

crea.lign Opaquer, 4 g



1 x GUM

crea.lign paste GUM, 3 g



1 x bright

1 x PL light

1 x PO opaque


1 x orange

1 x brown

**REF CLIGNSETN**

**crea.lign Start set N**


Single-hand grip crea.lign gel  
with UV LED



3 x

crea.lign gel, 5 g

10 x



crea.lign gel  
shades of your choice  
→ see page 16/17




REF CLIGNSET3A

crea.lign test set colours A2, A3, A3,5


Application cannula

1 x 10 

crea.lign paste dentine, 3 g

 1 x A2  
1 x A3  
1 x A3,5


Rotary piston

  
3 x






Pin cover

1 x 10 


crea.lign gel Enamel and Incisal,, 5 g

 1 x E2  
1 x E3  
1 x I1 opal  
1 x I2 blue


crea.lign Composite shaper kit

 1 x angle chisel  
 1 x cup chisel  
 1 x taper point  
 1 x cup round  
 1 x flat chisel

modelling liquid  
10 ml

1 x 


crea.lign gel dentine, 5 g

 1 x A2  
1 x A3  
1 x A3,5


visio.link

10 ml  
1 x 

crea.lign gel Transpa clear, 5 g

 1 x T1

Single-hand grip crea.lign gel  
with UV LED

1 x 


crea.lign gel GUM, 5 g

 1 x G1 light  
1 x G2 pink  
1 x G3 pink  
1 x G4 red

Unique brush no. 4

1 x 


Mixing pad

1 x 

MKZ primer  
5 ml

1 x 


crea.lign opaque, 4 g

 1 x A2  
1 x A3  
1 x A3,5


REF VLTOOLKIT

visio.lign tool kit finishing set for composite veneers

Holder for rotating tools and  
Abraso-Star gloss jar (glass jar is not filled!)

1 x 

diamond disc  
on both sides mounted,  
flexible, Ø 19 x 0.23 mm

1 x 

Ceragum medium lens  
Ø 22 x 4 mm

12 x 


Abraso-Fix round brush  
extra fine, yellow  
Ø 22 mm

8 x 


Hexagonal brush  
goat-hair white  
Ø 13 mm

15 x 

Cotton buff  
Ø 22 mm

15 x 


tungsten carbide burs  
Ø 4.0 mm – garnet

1 x 

tungsten carbide burs  
Ø 2.3 mm – pear

1 x 

crea.lign  
surface cleaner  
125 ml

1 x 


Glass jar  
unfilled for Abraso-Star glass

1 x 


tungsten carbide burs  
Ø 4.0 mm – garnet

1 x 


tungsten carbide burs  
Ø 2.3 mm – conical, round

1 x 

Quick-Mandrell,  
Gr. 2

1 x 

Abraso-Star glass  
high lustre polishing  
50 ml

2 x 

Acrypol  
polishing paste for veneering materials  
170 g

1 x 

3

---

**visio.paint**  
colour natural beauty

---



WITH THE VISIO.LIGN COLOR  
COMPOSITE PAINTS THE  
CHROMA OF A VENEER CAN  
BE OPTIMALLY ADAPTED WHILE  
REACHING PARTICULAR COLOR  
EFFECTS. WITH THE VISIO.LIGN  
SHIELD GLOSSY GLAZING  
AN OPTIMAL SURFACE GLOSS  
AND PROTECTION CAN BE  
OBTAINED.

STAINS  
GLOSSY GLAZING

# visio.lign | color



2,6 ml

REF  
**VLCOL**

## visio.lign color

visio.lign color are light-curing composite paints that can be used to effortlessly achieve colour effects on the surface of a veneer. The visio.lign color colours Body A-D can be used to optimally adjust a veneer's chroma. Furthermore, there is a wide range of visio.lign color paints available for efficient surface characterisation.

# visio.lign | shield

REF  
**VLSHIELDLV**



**LV** - Low Viscosity  
(liquid) 5 ml

REF  
**VLSHIELDHV**



**HV** - High Viscosity  
(viscous) 5 ml

## visio.lign shield

The new visio.lign shield gloss varnish for composite, PMMA and BioHPP allows you to create optimal surface gloss and protection. A surface finished with visio.lign shield gloss

varnish acts as an abrasion-resistant protective shield that won't discolour. Available in two viscosity variants - liquid LV and viscous HV.





**visio.lign** color Set

# visio.lign | sets

**REF VLCOLSET visio.lign color Set**

aluminium case,  
width 34 cm x height 25 cm x depth 4 cm



1 x

visio.lign shield LV  
5 ml

1 x



visio.lign shield HV  
5 ml

1 x



crea.lign  
surface cleaner  
30 ml

1 x

visio.lign color  
2,6 ml

je 1 x



Mixing pad

1 x



The assortment is composed of the following visio.lign colors:

- Body A
- Body B
- Body C
- Body D
- Grey-Blue
- Blue
- Yellow
- Brown
- Olive
- White
- Orange
- Black
- Rose
- Creme

Disposable brush

10 x



Brush holder

1 x



**REF VLCOLSKIT visio.lign color Starter Kit**

aluminium case,  
width 34 cm x height 25 cm x depth 4 cm



1 x

visio.lign shield LV  
5 ml

1 x



visio.lign shield HV  
5 ml

1 x



crea.lign  
surface cleaner  
30 ml

1 x



visio.lign color  
2,6 ml

je 1 x



Mixing pad

1 x



The assortment is composed of the following visio.lign colors:

- Body A
- Body B
- Grey-Blue
- Olive
- Orange
- Creme


Disposable brush

10 x



Brush holder

1 x



4





**novo.lign**  
veneers natural beauty

---

NOVO.LIGN VENEERS  
IMPRESS WITH THEIR UNIQUE  
PROPERTIES. WITH THE PRIMERS  
& BONDERS OF THE GROUP,  
THEY CAN BE ATTACHED TO  
ANY FRAMEWORK MATERIAL  
AND SUPPLEMENTED WITH  
CREA.LIGN. BECAUSE THEY  
HAVE THE SAME SHADES  
AND SHAPES, THEY CAN BE  
SEAMLESSLY COMBINED WITH  
NEO.LIGN DENTURE TEETH.

VENEERS

# novo.lign | Veneers

## FEATURES



The novo.lign veneers consist of high impact polymer composite filled with microceramic, and can be used for the permanent veneering of metal, ceramic and polymer frameworks. They are moulded from real rows of teeth and, therefore, give a perfect natural look. The morphological layer structure also ensures the veneer has a natural depth. Thanks to its exceptional mechanical properties,

such as its high flexural strength (140 MPa) and low modulus of elasticity (approx. 3,000 MPa) the physiological veneers have a shock-absorbing effect and are therefore, perfect for implant prosthetics. The homogeneous and sealed surface ensures long-term colour stability and plaque resistance. Because they have the same shades and shapes, they can be seamlessly combined with neo.lign denture teeth.

novo.lign veneers are available in anterior and posterior shapes.

Posterior shapes allow implementation of any occlusion concept through the multifunctional chewing surface design of G, W and L shapes



- The "extra" wide W designs are more suited to veneering bulky crowns and bridge frameworks than the "slimline" G designs...
- ...the L designs with "lingualised occlusal design" for manufacturing implant-supported full dentures.

All novo.lign anterior and posterior denture teeth are available in all classic A-D shades and BL3 bleached shade.

**novo.lign**

# novo.lign | Coding

All tooth cards/slides for the novo.lign veneers are printed with silver lettering.

## novo.lign A




**A** stands for anterior, which includes all anterior tooth designs for both the upper and lower jaws.

Tooth mould, e.g. D49

Tooth shade, e.g. A3

## novo.lign P

G designs




**P** stands for posterior, which includes all posterior tooth designs for both upper and lower jaws.

For novo.lign P variants of veneers, a distinction is made between G, W and L designs.

G, W and L designs are available in several sizes. Each size is assigned a numeric code. The size increases with the numeric value. The G designs of the veneers are available in sizes 3 and 4; the W designs in sizes 3, 4 and 5; the L designs in sizes 2 and 3.


W designs



➔ For sizing charts see page 37

G and W designs are grouped into quadrants, i.e. on a single tooth card/slide, for example, the lower jaw premolars and molars of the 3rd quadrant are all one design and size, e.g. W4.

L designs



On the tooth cards/slides for G, W and L designs, all four quadrants are shown in the form of four small boxes.

The corresponding quadrant is marked with a dot.

The L designs for the novo.lign P veneers are available in two sizes, L2 and L3. A tooth card/slide is fitted with either 6 upper jaw (16, 15, 14, 24, 25, 26) or 6 lower jaw posterior teeth (46, 45, 44, 34, 35, 36)

# novo.lign | Colours and forms


## novo.lign A

Anterior / front teeth



**Buccal wall thickness**  
**1.0 mm**

- 12 upper moulds
- 3 lower moulds



**Variety with system**

## novo.lign P

Posterior teeth

### G designs



**Buccal wall thickness**  
**1.2 mm**

- Multi-functional occlusal surface design for all occlusion concepts

Available in sizes




### L designs



**Buccal wall thickness**  
**1.0 mm**

- Lingualised occlusal design for full dentures

Available in sizes



### W designs




**Buccal wall thickness**  
**0.8 mm**

- Multi-functional occlusal design for all occlusion concepts
- Extra-wide tooth design, particularly recommended for veneering crowns and bridges

Available in sizes



All novo.lign A and novo.lign P shapes are available in the following colours:



# novo.lign | Range of designs

Upper anterior



C43

REF VOC43 \_\_\_ S



H46

REF VOH46 \_\_\_ S



A44

REF VOA44 \_\_\_ S



S46

REF VOS46 \_\_\_ S



F44

REF VOF44 \_\_\_ S



I47

REF VOI47 \_\_\_ S



I45

REF VOI45 \_\_\_ S



D48

REF VOD48 \_\_\_ S

Lower anterior



T35

REF VUT35 \_\_\_ S



D38

REF VUD38 \_\_\_ S

The colour completes the REF:



B L 3



A 1 0



A 2 0



A 3 0



A 3 5



A 4 0



B 1 0



B 2 0



M48

REF VOM48 \_\_\_ S



D49

REF VOD49 \_\_\_ S



B51

REF VOB51 \_\_\_ S



K53

REF VOK43 \_\_\_ S



D41

REF VUD41 \_\_\_ S

Upper/lower posterior



L2

REF VOWL2 \_\_\_




L2

REF VUWL2 \_\_\_



L3

REF VOWL3 \_\_\_



L3

REF VUWL3 \_\_\_

➔ For more posterior designs see page 38 and 39



B30



B40



C10



C20



C30



C40



D20



D30



D40

# novo.lign | Range of designs

Upper/lower posterior



W3

REF VO1W3 \_ \_ \_



W3

REF VO2W3 \_ \_ \_



W4

REF VO1W4 \_ \_ \_



W3

REF VU4W3 \_ \_ \_



W3

REF VU3W3 \_ \_ \_



W4

REF VU4W4 \_ \_ \_



G3

REF VO1G3 \_ \_ \_



G3

REF VO2G3 \_ \_ \_



G4

REF VO1G4 \_ \_ \_



G3

REF VU4G3 \_ \_ \_



G3

REF VU3G3 \_ \_ \_



G4

REF VU4G4 \_ \_ \_

The colour completes the REF:



B L 3



A 1 0



A 2 0



A 3 0



A 3 5



A 4 0



B 1 0



B 2 0





W4

REF VO2W4 \_ \_ \_



W5

REF VO1W5 \_ \_ \_



W5

REF VO2W5 \_ \_ \_



W4

REF VU3W4 \_ \_ \_



W5

REF VU4W5 \_ \_ \_



W5

REF VU3W5 \_ \_ \_



G4

REF VO2G4 \_ \_ \_



G4

REF VU3G4 \_ \_ \_

→ For more posterior designs see page 39



B30



B40



C10



C20



C30



C40



D20



D30




D40



# novo.lign | sets

## REF VLIGNSET2 visio.lign veneering system novo.lign A + P starter kit 2

visio.lign tooth cabinet / box




1 x

novo.lign A veneers - Anterior teeth




10 tooth cards/slides

novo.lign P veneers - Posterior teeth



8 tooth cards/slides

visio.link, 10 ml



1 x

Mixing cannula combo.lign, size 2



1 x 10

combo.lign luting composite, 8 g



1 x

Assortment in one shade of your choice

➔ Classic A to D shades and BL3 bleached shade

visio.lign tooth cabinet

Free selection of the color design

➔ see page 82

novo.lign A and P tooth cards/slides

Choice of anterior and posterior veneer designs

➔ see page 38 to 41


novo.lign P tooth cards/slides

Choice of posterior veneer designs

➔ For sizing charts see page 37

## REF VLIGNSET3 visio.lign veneering system starter kit 3

visio.lign tooth cabinet / box




1 x

novo.lign A veneers - Anterior teeth




13 tooth cards/slides

novo.lign P veneers - Posterior teeth




8 tooth cards/slides

Single-hand grip crea.lign gel with UV LED




2 x

visio.link, 10 ml



1 x

crea.lign opaque, 4 g



1 x

Pin cover




1 x 10

Application cannula



1 x 10

crea.lign gel, 5 g



2 x

Mixing cannula combo.lign, size 2



1 x 10

combo.lign luting composite, 8 g



1 x

Assortment in one shade of your choice

➔ Classic A to D shades and BL3 bleached shade

visio.lign tooth cabinet

Free selection of the color design

➔ see page 82

novo.lign A und P tooth cards/slides

Choice of anterior and posterior veneer designs

➔ see page 38 to 41

novo.lign P tooth cards/slides


Choice of posterior veneer designs

➔ For sizing charts see page 37

REF VLIGNSET4

## visio.lign Starterkit 4

visio.lign tooth cabinet / box



1 x

novo.lign A veneers - Anterior teeth



26 tooth cards/slides

novo.lign P veneers - Posterior teeth



24 tooth cards/slides

Single-hand grip crea.lign gel with UV LED



3 x

visio.link, 10 ml



1 x

Assortment in two shades of your choice

→ Classic A to D shades and BL3 bleached shade

visio.lign tooth cabinet

Free selection of the color design

→ see page 83

novo.lign A and P tooth cards/slides

Choice of anterior and posterior veneer designs

→ see page 38 bis 41

novo.lign P tooth cards/slides

Choice of posterior veneer designs


→ For sizing charts see page 37

crea.lign incisal gel and

choice of dentine shades

→ see page 16

crea.lign opaque, 4 g



2 x

Pin cover



1 x 10

Application cannula




1 x 10

crea.lign gel, 5 g




10 x

Mixing cannula combo.lign, size 2




1 x 10

combo.lign luting composite, 8 g




2 x

Mixing cannula green, size 4



1 x 12

visio.sil silicone for matrices, 50 ml



1 x

MKZ primer, 5 ml



1 x

modelling liquid, 10 ml



1 x


**REF VLIGNPRVN** Living mould chart **novo.lign** - novo.lign A and novo.lign P

visio.lign folding box



1 x

novo.lign A veneers - Anterior teeth




12 x

3 x

= 15 tooth cards/slides:  
1 tooth card/slide in each design

novo.lign P veneers - Posterior teeth




5 x

5 x

= 20 tooth cards/slides:  
each G and W design in sizes:  
G3, G4 and W3, W4, W5.

5 x

5 x



2 x

= 4 tooth cards/slides:  
each 1 x L design  
in sizes:  
L2 and L3


2 x

Assortment in one shade of your choice

➔ Classic A to D shades and  
BL3 bleached shade


**REF VLIGNPRVN** **novo.lign** shape kit - reduced quantity of shapes

visio.lign folding box



1 x

novo.lign A veneers - Anterior teeth




11 x

3 x

= 14 tooth cards/slides:  
1 tooth card/slide in each design  
(Upper jaw design K53 is not included)

novo.lign P veneers - Posterior teeth



5 x

5 x

= 20 tooth cards/slides:  
G and W design  
in sizes:  
G3 and G4,  
W3, W4 and W5.

5 x

5 x

Assortment in one shade of your choice

➔ Classic A to D shades and  
BL3 bleached shade

5



**neo.lign**

choose natural beauty

---

THE NEO.LIGN FULL TEETH **FULL TEETH**  
ARE IDENTICAL IN COLOUR,  
SHAPE AND MATERIAL TO  
THE NOVO.LIGN VENEERS.  
THIS MAKES THEM IDEALLY  
SUITED FOR AESTHETICALLY  
UNIFORM COMBINATION  
TECHNOLOGY. NEO.LIGN'S  
LINGUALISED CHEWING  
SURFACE DESIGN IS OPTIMAL  
FOR TOTALPROTHETICS.  
THE HIGH MECHANICAL  
STRENGTH OF THE  
NEO.LIGN TEETH MAKES  
THEM PARTICULARLY SUITABLE  
FOR IMPLANT-SUPPORTED  
DENTURES.

# neo.lign | Teeth



Natural mamelon structure with 4 chromatic layers. Additionally incorporated Transpa layer creates a "pearl effect" and a natural effect of depth.



Specially incorporated grooves reproduce the natural surface structure.



Optimised functional guide surfaces: The upper and lower jaw teeth interlock in a true-to-life manner.

Optimum curvature of the tooth. This supports the lip and produces a natural, youthful appearance.

The neo.lign full teeth consist of 5 thin layers of high-impact PMMA composite (polymethyl methacrylate) and 4 chromatic layers. They are suitable for the fabrication of full prosthetic restorations and for the combination technique. Due to their high mechanical strength, neo.lign teeth are perfectly suited for implant-supported dentures. Just like the novo.lign veneers, the neo.lign full teeth are based on natural teeth. This results in a particularly natural appearance. The morphological layer structure provides an additional and natural effect of depth.

The neo.lign full teeth feature the same shade, design and material as the novo.lign veneers. The layer structure of neo.lign and novo.lign is also perfectly coordinated. As a result, neo.lign and novo.lign can be combined without any problems.

A homogeneous and dense surface of the neo.lign full teeth ensures permanent shade stability and plaque resistance.

neo.lign full teeth are available as anterior and posterior teeth.

The posterior teeth enable the implementation of any occlusion concept thanks to the special occlusal designs of the neo.lign G and L designs.





- With the L designs, the visio.lign system offers a lingualised occlusal surface design especially for the fabrication of full prosthetic restorations.
- With the neo.lign G designs, a multifunctional occlusal surface design is available. Due to its slender shape, the G design is particularly suited for cases with limited space.

All neo.lign anterior and posterior designs are available in the classic A-D shades and in a bleach shade BL3.

## 4 CHROMATIC LAYERS

neo.lign

## neo.lign | Coding

All tooth cards/slides for the neo.lign teeth are printed with **gold lettering**.

### neo.lign A




**A** stands for anterior, which includes all anterior tooth designs for both the upper and lower jaws.

Tooth mould, e.g. **D48**

Tooth shade, e.g. **B4**

### neo.lign P

G designs



**P** stands for posterior, which includes all posterior tooth designs for both the upper and lower jaws.


For neo.lign P variants of posterior teeth, a distinction is made between **G and L designs**.

Both the G and L designs are available in **several sizes**. Each size is assigned a numeric code. The size increases with numeric value. G and L designs of full teeth are available in sizes 2, 3 and 4.

→ For sizing charts see page 51

The G shapes are grouped by quadrant, i.e. premolars and molars, e.g. of the lower jaw of the 4th quadrants of a shape and size, e.g. **G2**. On the toothed strips of the G and L shapes, all four quadrants are shown in the form of four small boxes. The corresponding quadrant is marked with a dot.

L designs



With the L-shapes, the toothed strips are equipped with either all upper jaw (16, 15, 14, 24, 25, 26) or all lower jaw teeth (46, 45, 44, 34, 35, 36).

# neo.lign | Colours and forms

## neo.lign A

Anterior teeth



- 13 upper moulds
- 4 lower moulds

Variety with system



## neo.lign P

Posterior teeth

### G designs



- Multi-functional occlusal design for all occlusion concepts

### L designs



- lingualised occlusal design for full dentures



Available in sizes




small  large

Available in sizes



small  large

All neo.lign A and neo.lign P shapes are available in the following colours:



# neo.lign | Variety of shapes

Upper anterior



F41

REF TOF41 \_\_\_ S



F44

REF TOF44 \_\_\_ S



B43

REF TOB43 \_\_\_ S



I45

REF TOI45 \_\_\_ S



C43

REF TOC43 \_\_\_ S



H46

REF TOH46 \_\_\_ S



A44

REF TOA44 \_\_\_ S



S46

REF TOS46 \_\_\_ S

Lower anterior



T33

REF TUT33 \_\_\_ S



T35

REF TUT35 \_\_\_ S

The colour completes the REF:



B 1 3



A 1 0



A 2 0



A 3 0



A 3 5



A 4 0



B 1 0



B 2 0



I47

REF TOI47 \_\_\_ S



B51

REF TOB51 \_\_\_ S



D48

REF TOD48 \_\_\_ S



M48

REF TOM48 \_\_\_ S



D49

REF TOD49 \_\_\_ S



D38

REF TUD38 \_\_\_ S



D41

REF TUD41 \_\_\_ S



B 3 0



B 4 0



C 1 0



C 2 0



C 3 0



C 4 0



D 2 0



D 3 0



D 4 0

# neo.lign | Variety of shapes

Upper/lower posterior



G2

REF TO1G2 \_ \_ \_



G2

REF TO2G2 \_ \_ \_



G3

REF TO1G3 \_ \_ \_



G2

REF TU4G2 \_ \_ \_




G2

REF TU3G2 \_ \_ \_



G3

REF TU4G3 \_ \_ \_




L1

REF TOWL1 \_ \_ \_



L2

REF TOWL2 \_ \_ \_



L1

REF TUWL1 \_ \_ \_



L2

REF TUWL2 \_ \_ \_

The colour completes the REF:



B 1 3



A 1 0



A 2 0



A 3 0



A 3 5



A 4 0



B 1 0



B 2 0



G3

REF TO2G3 \_ \_ \_



G4

REF TO1G4 \_ \_ \_



G4

REF TO2G4 \_ \_ \_



G3

REF TU3G3 \_ \_ \_




G4

REF TU4G4 \_ \_ \_




G4

REF TU3G4 \_ \_ \_




L3

REF TOWL3 \_ \_ \_




L4

REF TOWL4 \_ \_ \_



L3

REF TUWL3 \_ \_ \_



L4

REF TUWL4 \_ \_ \_



B30



B40



C10



C20



C30



C40



D20



D30



D40






# neo.lign | sets


**REF V LIGNPRVT** Living mould chart **neo.lign** - neo.lign A and neo.lign P

visio.lign folding box



1 x

neo.lign A full teeth - Anterior teeth




13 x

4 x

= 17 tooth cards/slides:  
1 tooth card/slide in each design

neo.lign P full teeth - Posterior teeth G designs



3 x


3 x

3 x

3 x

= 12 tooth cards/slides:  
each 1 x G design in sizes:  
G2, G3 and G4

neo.lign P full teeth - Posterior teeth L designs



4 x

4 x


= 8 tooth cards/slides:  
each 1 x L design in sizes:  
L2, L3 and L4

Range of all tooth shapes in a colour of your choice

➔ Classic A to D shades and BL3 bleached shade

**REF V LIGNPRVT** **neo.lign** shape kit - reduced quantity of shapes

neo.lign A full teeth - Anterior teeth




13 x

4 x

= 17 tooth cards/slides:  
1 tooth card/slide in each design

neo.lign P full teeth - Posterior teeth L designs



4 x

4 x

= 8 tooth cards/slides:  
each 1 x L design in sizes:  
L2, L3 and L4

Assortment in one shade of your choice

➔ Classic A to D shades and BL3 bleached shade


6

---

# **bond.lign**

secure natural beauty

---



BOND.LIGN IS A UNIQUE COMPOSITE SYSTEM FOR PERMANENTLY SECURING THE ADHESION OF ALL VISIO.LIGN SYSTEM COMPONENTS AND COMPOSITES TO ALL COMMON FRAMEWORK MATERIALS SUCH AS NEM, EM, TITANIUM, HIPC (POLYMERS/ COMPOSITES), PEEK (E.G. BIOHPP) AND ZIRCONIUM OXIDE.

PRIMER AND BONDER  
OPAQUE  
LUTING COMPOSITE

## Primer and bonder



REF **MKZ02004**

5 ml

### MKZ primer

Adhesive bond of composite to metal frameworks

- Titanium
- Zirconium oxide
- Non-precious metal alloys
- Oxide ceramics

The MKZ primer is used to condition all metallic NEM alloys, titanium, zirconium oxide and oxide ceramics in order to create a chemical bond to composites such as crea.lign or novo.lign as well as to all other plastics (e.g. top.lign, PMMA, etc.). The MKZ primer can be used intraorally and is ideally suited for bonding of individual abutments.



REF **MKZEM004**

4 ml

### MKZ EM activator

Adhesive bond of composite to precious metal frameworks

- Precious metal alloys (Au, Ag, Pt, Pd)
- eco-alloys  
(low precious metal alloys)

The MKZ EM activator in combination with the MKZ primer is used to condition precious metal frameworks (Au, Ag, Pt, Pd), in order to create a chemical bond to composites such as crea.lign or novo.lign. Like the MKZ primer, the MKZ EM activator can be used intraorally.

**⚠** Only use in combination with MKZ primer!



#### Application MKZ EM activator

1:1 mix with MKZ primer.

REF **APK25003**

4 ml

## K primer

Adhesive bonding of composite to oxide ceramic frameworks

- Lithium (di)silicate
- Veneering and press ceramics

The K primer produces the adhesive bond between lithium (-di) silicate, veneering ceramics and press ceramics to composites such as crea.lign or novo.lign. Ideal when repairing ceramic veneers, as the K primer may be used intraorally.



### silanisation

K primer is also suitable for silanisation of surfaces.

REF **VLPMA10**

10 ml

## visio.link

Adhesive bonding of composite to plastics

- PMMA denture resins
- Composite  
(veneer composites/composite teeth)
- High-performance polymers BioHPP
- PEEK/PEKK/PAEK

The light-curing visio.link is used to create the bond between PMMA prosthesis plastics, high-impact PMMA composite materials, high-performance polymers such as BioHPP and materials such as PEEK/PEKK/PAEK to other composites such as crea.lign or novo.lign. visio.link can also be used intraorally.

# bond.lign

## Opaquer combo.lign

light



OPL

A1-A3

B1-B2

C1-C2

medium



OPM

A3.5

B3-B4

D2-D3

intensive



OPI

A4

C3-C4

D4

GUM



OPG

Gingiva



REF CO1X4

4 g

### Opaque combo.lign

is a light-curing and self-curing opaque that guarantees safe curing. It is colour-coordinated with the use of combo.lign luting composite and novo.lign veneers. The university-tested composite system is compatible with all metal primers and silanisation processes and is ideal when using mechanical retentions. The combo.lign opaque is available in one gingival and three tooth colours. In combination with combo.lign luting composite, this guarantees colour stability in all A-D colours.



### combo.lign opaque catalyst

Light polymerisation is necessary to achieve an optimal bond!

For self-hardening, use a 1:1 mix of combo.lign opaque catalyst and combo.lign opaque mass, in four different shades.



REF CO1X4KAT


4 g

### Opaquer combo.lign Catalyst

### Shade combination table combo.lign opaque

combo.lign opaque				
system shades	light	medium	intensive	GUM
shades	A1 - A3 / B1 - B2 / C1 - C2	A3.5 / B3 - B4 / D2 - D3	A4 / C3 - C4 / D4	Gingiva shades

# crea.lign opaque



**⚠** With pre-coloured zirconium frameworks, the colour shines through, which means that conventional opaque does not have to be applied - provided the colour of the framework corresponds to the desired (A-D) colour.



REF CLFHO 4 g

### crea.lign opaque

The light-curing crea.lign opaque provides shade reliability in the freestyling layer and when using novo.lign veneers. Together with the MKZ primer, it forms the basis for optimum bonding with metallic and ceramic framework materials and therefore for long-lasting prosthetic restorations without marginal gaps. All shades of the classic A-D system are covered by just nine opaques. The special GUM opaque is suitable for all GUM shades.



transparent

REF CLFHOPZ 4 g

### crea.lign opaque Z, zirconium liner


The light-curing zirconium liner is a transparent opaque. When using zirconium it is essential to apply zirconium liner after applying the appropriate primer (MKZ primer) to achieve reliable bonding of zirconium oxide and composite.

## Shade combination table crea.lign opaque

crea.lign opaque system shades	1	2	3	4	5	6	7	8	9	GUM
shades	A1 / B2	A2	A3	B1 / C1 / BL3	C2 / C3 / D2 / D4	B3 / B4	A3.5	A4 / C4	D3	Gingiva shades

The transparent crea.lign® opaque is suitable for pre-coloured frameworks, e.g. zirconium.


# combo.lign luting composite



REF **CO2X4** \_ \_ \_  
8 g

combo.lign is a light-curing and dual-cure luting composite for reliable colour reproduction and optimal bonding between novo.lign veneer and all framework materials.

# DTK adhesive



The DTK dual-cure adhesive enables a durable and secure bond between all prosthetic materials. Ideal for bonded, individual abutments.

Thanks to its sterilisability, contamination is avoided and maximum hygiene is possible. The DTK adhesive is suitable for intraoral and extraoral use.





**bond.lign**  
full range bonding kit




**bond.lign**  
chipping repair kit

# bond.lign | Sets

REF CLFHOPSET

## crea.lign Opaquer set

crea.lign opaque, 4 g




crea.lign opaque

→ see page 63

REF COMBOKIT

## combo.lign complete kit

combo.lign, 8g



→ see page 64


REF 54001285

## DTK adhesive opak assortment (5 Parts)

Mixing cannulas



1 x 10

Double-mix cartridge  
DTK adhesive opaque, 8 g

1 x

Syringe plunger

1 x

Disposable brush



1 x 10

Brush holder



1 x

DTK adhesive

→ see page 64


REF 54001185

## DTK adhesive assortment (5 parts)

Mixing cannulas



1 x 10

Double-mix cartridge  
DTK adhesive transparent, 8 g

1 x

Syringe plunger

1 x

Disposable brush



1 x 10

Brush holder



1 x

DTK adhesive

→ see page 64

REF 54001295


## refill pack

Double-mix cartridge  
DTK adhesive opaque, 8 g

1 x

REF 54001195

## refill pack


Double-mix cartridge  
DTK adhesive transparent, 8 g

1 x

**REF 5400bond** Full range bonding kit

bredent medical:  
**REF 5800bond**

Plastic case




Mixing pad



MKZ primer




K primer




visio.link



Silano-Pen plastic cups



One-time dental brush



Brush holder, black




Intra-oral-tips + Mixing cannulas



DTK adhesive opaque, 8 g



DTK adhesive transparent, 8 g



FGP insulation, 3 ml



Fixation screw M1.4



**REF CLIGNSETA** Chipping repair kit

Aluminium case



MKZ primer



MKZ EM activator



K primer



crea.lign paste dentine, 3 g



Rotary piston



Silano-Pen plastic cups




Pin cover



Application cannula



crea.lign gel incisal, 5 g




crea.lign opaque, 4 g




visio.link




crea.lign gel dentine 5 g



crea.lign opaque Z, zirconium liner 4 g



Processing instructions insert



dentist

7



**visio.CAM**














design natural beauty

---

THE MATERIAL COMPETENCE  
OF THE VISIO.LIGN SYSTEM IS  
ALSO AVAILABLE FOR  
CAD/CAM PROCESSING.  
THE BRECAM.HIPC BLANKS  
OF THE VISIO.CAM LINE ARE  
COMPATIBLE WITH THE ENTIRE  
VISIO.LIGN SYSTEM.

COMPOSITE BLANKS

## Composite blanks

breCAM.HIPC	W DCS		ZZ
	Ø 98,5 mm		Ø 95 mm
	16 mm	20 mm	20 mm
 <b>bleach</b>	REF 54003418	REF 54003398	REF 54003438
 <b>A1</b>	REF 54003419	REF 54003399	REF 54003439
 <b>A2</b>	REF 54003420	REF 54003400	REF 54003440
 <b>A3</b>	REF 54003421	REF 54003401	REF 54003441
 <b>A3,5</b>	REF 54003422	REF 54003402	REF 54003442
 <b>A4</b>	REF 54003423	REF 54003403	REF 54003443
 <b>B2</b>	REF 54003424	REF 54003404	REF 54003444
 <b>B3</b>	REF 54003425	REF 54003405	REF 54003445
 <b>C2</b>	REF 54003426	REF 54003406	REF 54003446
 <b>C3</b>	REF 54003427	REF 54003407	REF 54003447
 <b>D2</b>	REF 54003428	REF 54003408	REF 54003448
 <b>D3</b>	REF 54003429	REF 54003409	REF 54003449
 <b>transparent</b>	REF 54003430	REF 54003410	REF 54003450

The advantages of the visio.lign system have also been transferred to digital processing with visio.CAM. HIPC means "high impact polymer composite" and originates from the development of the visio.lign system. Chemically, breCAM.HIPC corresponds to novo.lign veneers, making this material ideal for permanent dentures. The production takes place with pressure and heat, similar to prefabricated teeth made of plastic, at approx. 120 degrees and 250 bar pressure.



The absence of dental glass and light-curing plastic ensures high colour stability and plaque resistance. The breCAM.HIPC composite blanks of the visio.CAM line are optimally matched to the system components of the visio.lign system.

They are also available in 20 mm height for DCS, Wieland and zirconium tooth machines and 16 mm height for DCS and Wieland machines. Available in the most common A-D colours.

**visio.CAM**

8



# Tools



support natural beauty

---

THE FOLLOWING CHAPTER PRESENTS DEVICES, TOOLS AND ACCESSORIES THAT ARE RECOMMENDED FOR OPTIMAL AND EASY USE OF VISIO.LIGN PRODUCTS.

DEVICES,  
TOOLS  
AND ACCESSORIES

# Tools | bre.Lux

bre.Lux LightPen  
Set  
Angled hand light  
+ Anti-glare shield (orange)  
+ Adapter  
+ Charging cradle

REF 1400204

1 Set



## Technische Daten bre.lux Light Pen

Light source	blue light
Light range	385 - 515 nm
Light intensity	1000 mW / cm <sup>2</sup> - 1700 mW / cm <sup>2</sup>
Rechargeable lithium batteries	Model ICR 18650
Battery voltage and capacity	3,6 V / 2000 mAh
Adapter input	100 V - 240 V / 50Hz - 60 Hz

The device complies with the EMC\* Directive 2014/30/EU.

\*Electromagnetic compatibility

bre.Lux Power Unit 2  
light tray including 30 sticks

REF 14001002

1 x 30 sticks

bre.Lux Power Unit 2  
light tray sticks cylinder

REF 14001003

1 x 30 pieces

bre.Lux Power Unit 2  
light tray sticks taper and cone

REF 14001004

1 x 28 pieces

The LED light curing device  
bre.Lux Power Unit 2 set

impresses with safety and speed in the polymerisation of objects. The innovative LED technology uses the completely relevant light spectrum of 400 nm and thus achieves a high quality of polymerisation. Different wavelengths reach different depths - this enables consistent and faster curing. Freely programmable for any requirement, including new materials.

## Technical data

### bre.Lux PowerUnit 2

#### Basic unit

U	100 - 240 VAC
P	130 W
Frq	50/60 Hz
Fuse	T 2.0 A
Light range	370 - 500 nm

## Advantages

### Reliable polymerisation

- all light waves available  
= maximum penetration  
= top quality of the result

### Unsurpassed polymerisation speed

- 72 lights + full range  
+ rotary plate + light tray

### Reliability

- Careful and gradual increase of luminous power
- Temperature control 45°-55° C - up to 65° C possible for individually adjustable programmes
- Overheat protection
- Automatic deactivation of the light with acoustic signal at the end of the programme

### bre.Lux PowerUnit 2 Set Basic unit and hand lamp including accessories

- bre.Lux Power Unit 2 basic unit
- + bre.Lux LED N2 hand lamp  
(including flexholder, plug-in power unit)
- + Rotary plate
- + Light tray (including 30 sticks)

REF 14001015

1 piece

### bre.Lux Power Unit 2 basic unit including accessories Full range system

REF 14001000

1 piece



- 10 individually adjustable programmes +  
1 continuous light programme

#### Long service life of the LED lights

- 20,000 operating hours or 12 years

#### Ease of use

- Easy and fast access to programmes
- Programmable to satisfy individual requirements
- Spacious drawer
- Extremely silent and compact unit



→ polymerisation times bre.Lux Power Unit 2  
see pages 76/77

# Polymerisation times

## bre.Lux PowerUnit 2

Manufacturer	Product name	bre.Lux PowerUnit 2 curing times in seconds [s] / Luminous power in percent [%]				
		bre.Lux LightPen (Angled hand light)		bre.Lux PowerUnit2 (basic unit)		
		Partial polymerisation (fixation/ partial curing)	Intermediate polymerisation (intermediate curing of individual layers)	Partial polymerisation (fixation/ partial curing)	Intermediate polymerisation (intermediate curing of individual layers)	Final polymerisation (final end hardness of material)
bredent	visio.link	N/A	30 s	N/A	90 s	90 s
	combo.lign bonding composite	30 s	N/A	20 s / 50 % (iProg)	180 s	180 s
	crea.lign gel	15 s	15 s	20 s / 50 % (iProg)	180 s	360 s
	crea.lign paste	15 s	15 s	20 s / 50 % (iProg)	180 s	360 s
	combo.lign Opaker	N/A	N/A	N/A	180 s	360 s
	crea.lign Opaker	N/A	N/A	N/A	180 s	360 s
	visio.lign color/visio.lign shield	10 s	N/A	N/A	90 s	180 s
	crea.lign Stains	N/A	N/A	N/A	180 s	360 s
	Ropak UV	N/A	N/A	N/A	180 s **	360 s
	Ropak Compact opaquer UV	N/A	N/A	N/A	180 s **	360 s
	Ropak Compact opaquer tooth-coloured UV	N/A	N/A	N/A	180 s	360 s
	compoForm UV	15 s	N/A	20 s / 50 % (iProg)	90 s	180 s
	Tray material UV	N/A	N/A	40 s / 50 % (iProg)	90 s	180 s (1)
	Die varnish, light-curing	N/A	30 s *	20 s / 50 % (iProg)	90 s	180 s
	SERACOLL UV	N/A	15 s	N/A	20 s / 50 % (iProg)	40 s / 50 % (iProg)
Qu-connector	N/A	30 s	N/A	90 s	90 s	
Heraeus	Signum	N/A	N/A	20 s / 50 % (iProg)	180 s	360 s
	Palatray XL	N/A	N/A	40 s / 50 % (iProg)	90 s	360 s
Shofu	Solidex	N/A	N/A	20 s / 50 % (iProg)	180 s	360 s
GC	Gradia	15 s	N/A	20 s / 50 % (iProg)	180 s	360 s
Wegold	S-Lay	N/A	N/A	20 s / 50 % (iProg)	180 s	360 s
VITA	VITA VM LC PRE OPAQUE / VITA VM LC OPAQUE PASTE	N/A	N/A	N/A	N/A	180 s
	VITA VM LC OPAQUE Pulver	N/A	N/A	N/A	360 s	360 s
	VITA VM LC Composite	30 s	N/A	40 s / 50 % (iProg)	180 s	360 s (2)
Degudent	Degudent injoy	N/A	N/A	20 s / 50 % (iProg)	180 s	360 s

### Notes:

- \* In case of one order only.
- \*\* Only apply opaquer in two coats.
- (1) For UV tray material, upper and lower side each 1 x 180 s.
- (2) For pontics, layer thickness up to max. 2 mm.
- (iProg) Individual programming necessary: See new parameters, including the bonding of prefabricated wax parts. Please create an individual programme with 50% luminous power without step function (heat reduction!)
- N/A Not applicable.

The polymerisation times are guide values for intact devices.

**Individual programmes:** Metal frameworks store the thermal energy of light more than pure polymers. Heat can have a positive effect on materials in the form of post-treatment or, in the case of too much heat, can lead to embrittlement or stresses. The development of heat can be conveniently controlled by the user by automatically adapting the stored (default) programs to the desired material conditions. For metal-free restorations or constructions with material thicknesses of more than 2mm, the option up to 100% power "Red. Power off" is recommended. For restorations that include metal components or involve high material shrinkage, it is recommended to reduce the light power: "Red. Power on". However, raising the power to 100% is always possible in the individual programs without having to re-adjust the unit.  
In some cases, the curing times may change proportionally.

**The hand lamp can be used as an alternative to the bre.Lux PowerUnit 2 for partial or intermediate polymerisation. The final curing is always done using the bre.Lux PowerUnit 2!**

## Other suitable light polymerisation devices

Polymerisation times for the visio.lign system components visio.link, combo.lign and crea.lign

Manufacturer	Product name	Wavelength [nm]*	Polymerisation times in seconds [s]			
			visio.link	combo.lign	crea.lign / crea.lign - Opaker / combo.lign - Opaker	
bredent	bre.Lux PowerUnit 2	370 - 500 nm	90 s	180 s	360 s	
Dentsply	Triad 2000	400 - 500 nm	180 s	360 s	600 s	
Degudent	Eclipse	k.A.	60 s	180 s	360 s	
Heraeus Kulzer	Dentacolor XS, Uni XS, Heraflash	320 - 520 nm	90 s	180 s	360 s	
GC	Labolight LV-III	380 - 490 nm	120 s	300 s	600 s	
Ivoclar Vivadent	Targes Power Ofen, Luminat 100	400 - 580 nm	240 s	180 s	480 s	
Schütz Dental	Spektra 200	310 - 500 nm	120 s	180 s	360 s	
Shofu Dental	Solidilite	400 - 500 nm	90 s	180 s	360 s	
Kuraray Dental	CS 110	k.A.	120 s	300 s	480 s	
Hager & Werken	Speed Labolight	320 - 550 nm	90 s	180 s	480 s	
3M ESPE	Visio Beta	400 - 500 nm	neu: P1 - P4	> 240 s (P2)	420 s (P2)	900 s (P1)
			alt: U0 - U3	420 s (U1, U3)	900 s (U0)	900 s (U0)

### Notes:

\* Manufacturer's data  
n/a No information

## Layer thickness and shade combination tables

Polymerisation times for specific layer thicknesses

Material	max. layer thickness [mm]	bre.Lux PowerUnit 2 polymerisation times in seconds [s]	
		Intermediate polymerisation (intermediate curing of individual layers)	Final polymerisation (final end hardness of material)
crea.lign Enamel	1 mm	180 s	360 s
crea.lign Incisal	1 mm	180 s	360 s
crea.lign Transpa clear	1 mm	180 s	360 s
crea.lign Dentin	1 mm	180 s	360 s
crea.lign Opaque Dentin	1 mm	180 s	360 s
crea.lign High Translucent	1 mm	180 s	360 s
crea.lign Modifier	1 mm	180 s	360 s
crea.lign GUM	1 mm	180 s	360 s
crea.lign Stains	0,3 mm	180 s	360 s
crea.lign paste	2 mm	180 s	360 s
combo.lign	2 mm	180 s	180 s
crea.lign Opaker	0,1 mm	180 s	360 s
combo.lign Opaker	0,1 mm	180 s	180 s

# Tools

## visio.lign Toolkit

The combination of the materials of the visio.lign veneering system and the visio.lign toolkit results in a surface featuring plaque resistance and shade stability, and its quality and resistance is identical to that of a ceramic material.

### Composite surfaces like ceramics

The visio.lign toolkit has been designed specifically for the finishing of composites, such as visio.lign veneers, and provides tools and polishing pastes for a perfect finish.

Tool stand and glass jar without tools



REF **VL1SHELL**

1 piece, 2 parts

REF **VLTOOLKIT**

1 piece, 14 parts

Preparation set for composite veneers

## Benefits

### of the visio.lign toolkit

- Burs with relief produce smooth surfaces and reduce polishing time
- Prepolishing and high lustre polishing pastes are included in the set
- Tools are arranged at different heights to ensure easy access and grasp
- Removable glass pots to avoid drying out of high lustre polishing paste
- Pictograms and reference numbers are printed on the tools and stand for better visual control
- Three empty spaces for additional tools

### crea.lign surface cleaner

Cleaner for removing the dispersion layer



REF **43000600**

125 ml

Abraso-Star gloss high lustre polishing



REF **52000163**

2 x 50 ml

Acrypol polishing paste for veneering materials



REF **52000170**

170 g

# crea.lign gel

Stainless steel syringe holder  
for crea.lign syringe



Stainless steel syringe holder for up to 6 crea.lign syringes including the single-hand grip.

REF 32000960

1 piece

Syringe holder for  
crea.lign 5 g and stains



Up to 12 crea.lign syringes including the single-hand grip can be housed in this crea.lign syringe holder.

REF 32000443

1 piece

Application cannulas  
+ pin cover



Lightproof protective cover, suitable for all application cannulas

Light-blocking application cannulas suitable for all crea.lign 5 g syringes

REF 32000940

10 pieces each

Rotary piston  
with UV-LED



black with UV-LED

Single-hand grip for all crea.lign gel 5 g syringes, for simple and controlled dispensing

REF 32000445

1 piece

Mixing pad 100 x100 mm



REF 32000460

1 piece, 100 sheets

Mixing pads  
35 x 50 x10 mm



REF 33001144

10 pieces, 100 sheets each



Single-hand grip +  
application cannula  
for crea.lign gel syringes



The single-hand grip has been specially designed for crea.lign gel syringes and, in combination with the application cannulas, makes it easier for the user to ensure precise dosing with one hand. Using the single-hand grip, the syringe sits firmly in the palm of the hand and the liquid gel can be dispensed by turning the ribbed rotating device using the thumb or index finger of the same hand.

The other hand remains free to carry out other work, e.g. holding the object being processed.



The UV LED hand lamp is integrated into the single-hand grip and can be removed as required. The applied crea.lign gel can be quickly hardened with the aid of this.

# Tools

## crea.lign paste

syringe holder for  
crea.lign paste



Syringe holder in transparent plastic for five crea.lign pastes, unloaded.

REF 32000447

1 piece

crea.lign  
composite shaper kit



The modelling instruments made of silicone have been specially developed for the crea.lign paste. They serve as application and modelling aids.

REF 32000480

1 piece

Rotary piston  
crea.lign paste



Rotary piston for all crea.lign pastes, for simple and controlled dispensing

REF 32000444

1 piece

## novo.lign

Thermo-Pen



The Thermo-Pen is a flameless hot-air device with piezo technology. The veneers can be heated without the risk of burning them. The temperature required to shape the veneers is approx. 250 C. Expansion improves coating of telescopic crowns, for example.

REF 11001470

1 piece



We recommend the Thermo-Pen  
for thermoplastic shaping

1




Heat veneer

2



Bend hot veneer upwards





# bond.lign

bond.lign stainless steel tripod  
(unloaded)



High-quality metal stand with removable tubs - can be used on both sides (large and small cooling). Brush holder and spacious storage.

REF 32000970

1 piece

Dispenser for 5 ml 1:1  
combo.lign luting composite



⚠ Cartridge not included within the scope of delivery.

Dispenser for all combo.lign cartridges, for simple and controlled dispensing of combo.lign luting composite.

REF 32000441

1 piece

combo.lign  
mixing cannulas  
size 2



Lightproof mixing cannulas for all combo.lign luting composites

REF COMKG210

1 x 10 pieces

# Preparatory work

Multisil dosing device  
visio.sil silicone dispenser



⚠ Cartridge not included within the scope of delivery.

Dispenser for visio.sil, visio.sil ILT and visio.sil fix cartridges, for simple and controlled dispensing of the silicone.

REF 32000440

1 piece

Mixing cannulas size 2 yellow  
for visio.sil fix silicone



Mixing cannula for high-precision silicone for matrices.

REF 32000451

1 x 12 pieces

Mixing cannulas size 4 green  
for visio.sil and visio.sil ILT silicones



Mixing cannulas for transparent visio.sil silicone.

REF 32000457

1 x 12 pieces

# Tools

## visio.lign tooth cabinet / box

including 10 push inserts, unloaded



REF VLIGNBOX3

1 piece

Tooth and tray cabinet for storing neo.lign / novo.lign anterior and posterior sets

REF VLIGNBOX2

1 piece

### Loading drawer inserts for



26 x neo.lign/novo.lign P posterior sets (8/6 toothed strips) + 13 x neo.lign or novo.lign A/P sets

REF 00745961

1 piece



52 x neo.lign/novo.lign A anterior sets/52 x neo.lign/novo.lign P posterior sets (set of 4)

REF 00745962

1 piece



9 x crea.lign opaque and / or crea.lign paste syringe

REF 00745824

1 piece



6 x combo.lign opaque syringe

REF 00745827

1 piece



10 x crea.lign gel syringe  
3 x single-hand grip

REF 00745847

1 piece



8 x crea.lign gel and/or crea.lign paste syringe  
2 x adhesive bond liquid vials  
1 x silicone cartridge

REF 00745816

1 piece

# Material processing matrix

Material	Processing			
	Veneering	Freestyle layering	Freestyle layering with Gnathoflex	Freestyle layering in a matrix
visio.sil fix	•			
visio.sil	•			
visio.sil ILT	•			•
MKZ primer	•	•	•	•
MKZ EM activator	•	•	•	•
K primer	•	•	•	•
visio.link	•	•	•	•
bre.Lux Power Unit 2	•	•	•	•
Gnathoflex			•	
crea.lign modelling liquid	•	•	•	•
combo.lign opaque	•			
crea.lign opaque	•	•	•	•
combo.lign luting composite	•			
novo.lign Veneers	•			
crea.lign gel	•	•	•	•
crea.lign paste	•	•	•	•
visio.lign color/visio.lign shield	•	•	•	•
visio.lign Toolkit	•	•	•	•

# Media



## Preparatory work

REF 000482GB

Inverse Layering Technique  
*(by Vincenzo Musella)*



## visio.lign

REF 000224GB

brochure  
visio.lign

REF 000234GB

brochure  
visio.lign Instructions for use

REF 000274GB

visio.lign Binder  
english empty

REF 0091161GB

brochure  
visio.lign modern art

REF 0091177GB

brochure  
modern art Manual



## crea.lign

REF 000577GB

brochure  
crea.lign

REF 000651GB

brochure  
red-white aesthetic set

REF 009907GB

Flyer  
crea.lign  
paste

REF 009833GB

brochure  
crea.lign freestyle manual

REF 009485GB

Inlay  
Layering scheme for RED-  
WHITE individualization

REF 992990GB

Dental aesthetics book of  
Vincenzo Musella  
english

REF 009837GB

crea.lign  
Layering instructions



## novo.lign

REF 009939GB

brochure  
novo.lign

REF 000202GB

brochure  
range of designs

REF 000754GB

brochure  
Quick Guide SKY top.lign  
(& novo.lign)



neo.lign

REF 000305GB

brochure  
neo.lign

REF 000329GB

Design chart  
neo.lign Full teeth

REF 000468GB

brochure neo.lign  
Set-up instructions

REF 000226GB

neo.lign Curriculum  
Modern Prosthetics

REF 0091006GB

neo.lign  
Full teeth



bond.lign

REF 009539GB

Indication  
of the visio.lign primers

REF 009893GB

Flyer  
bond.lign/DTK adhesive

REF 000477GB

Data Sheet  
MKZ Primer / K-Primer

REF 009894GB

Adhesive protocol  
Full Range Bonding Kit



visio.paint

REF 009681GB

Quick reference card  
visio.lign color &  
visio.lign shield

REF 0091153GB

Material Report  
visio.lign shield &  
visio.lign color



visio.CAM

REF 000500GB

Material catalogue  
CAD/CAM



Tools

REF 000234GB

Prospekt  
visio.lign  
Manual

REF 000590GB

brochure  
bre.Lux PowerUnit 2

REF 0091174GB

brochure  
bre.Lux LightPen

# Index

## A

Abraso-Star high lustre polishing.....	78
Accessories	
Application cannulas + pin cover.....	79
bond.lign stainless steel tripod.....	81
combo.lign mixing cannulas.....	81
crea.lign composite shaper kit.....	80
Dispenser combo.lign luting composite.....	81
Mixing cannulas for visio.sil and visio.sil ILT silicones.....	81
Mixing cannulas for visio.sil fix silicone.....	81
Mixing pad.....	79
Rotary piston crea.lign paste.....	80
Stainless steel syringe holder for crea.lign syringe.....	79
Syringe holder for crea.lign 5 g & stains.....	79
syringe holder for crea.lign paste.....	80
Thermo-Pen.....	80
Tool stand and glass jar, without tools.....	78
visio.sil silicone dispenser Multisil dosing device.....	81
Acrypol polishing paste for veneering materials.....	78
adhesive DTK adhesive.....	64
Application cannulas + pin cover.....	79
visio.lign color.....	28
visio.lign shield.....	28

## B

Bonder visio.link.....	61
bre.Lux	
LightPen.....	74
Power Unit 2 basic unit.....	75
Power Unit 2 light tray including 30 sticks.....	74
Power Unit 2 light tray sticks cylinder.....	74
Power Unit 2 light tray sticks taper and cone.....	74
Power Unit 2 set.....	75
breCAM.HIPC composite blanks visio.CAM.....	70
Brush finishing set for composite veneers visio.lign toolkit.....	25, 78

## C

Chipping repair kit.....	67
cleaner crea.lign surface cleaner.....	78
combo.lign luting composite.....	64
complete kit.....	66
Mixing cannulas.....	81
Opaque catalyst.....	62
Opaque.....	62
Composite blanks visio.CAM breCAM.HIPC.....	70
crea.lign	
Assortment with syringe holder.....	22
Composite shaper kit.....	80
freestyle intro kit.....	23
gel dentine.....	16
gel GUM.....	16
gel Incisal.....	16
gel Modifier.....	16
gel Transpa Clear.....	16
gel.....	12, 16
Modelling liquid.....	18
Opaque set.....	66
Opaque Z, zirconium liner.....	63
Opaque.....	63
paste complete kit.....	23
paste dentine.....	17
paste GUM.....	17
paste Incisal.....	17
paste starter kit, 10 x crea.lign paste.....	22

paste starter kit, 9 pastes can be selected .....	22
paste Transpa opal .....	17
paste .....	12, 17
stainless steel syringe holder .....	79
Stains .....	16
Start set N .....	24
surface cleaner .....	78
Syringe holder for crea.lign 5 g and stains .....	79
syringe holder stainless steel .....	79
test set colours A2, A3, A3,5 .....	25
veneering system paste + gel, complete set .....	21
Cutter finishing set for composite veneers visio.lign toolkit .....	25, 78

## D

Dispenser	
combo.lign luting composite .....	81
visio.sil silicones Multisil dosing device .....	81
Dosing device	
Dispenser combo.lign luting composite .....	81
visio.sil silicone dispenser Multisil dosing device .....	81
DTK adhesive .....	64
abutment assortment .....	66
assortment transparent .....	66
opaque abutment set .....	66
opaque assortment .....	66

## F

finishing set for composite veneers visio.lign toolkit .....	25, 78
Full range bonding kit .....	67

## K

K primer .....	61
----------------	----

## L

Light curing device	
bre.Lux LightPen .....	74
bre.Lux Power Unit 2 basic unit .....	75
bre.Lux Power Unit 2 light tray including 30 sticks .....	74
bre.Lux Power Unit 2 light tray sticks cylinder .....	74
bre.Lux Power Unit 2 light tray sticks taper and cone .....	74
bre.Lux Power Unit 2 set .....	75
Liquid crea.lign modelling liquid .....	18
luting composite	
combo.lign luting composite .....	64
combo.lign mixing cannulas .....	75
Dispenser combo.lign luting composite .....	75

## M

Mixing cannulas	
combo.lign mixing cannulas .....	81
Mixing cannulas for visio.sil and visio.sil ILT silicones .....	81
Mixing cannulas for visio.sil fix silicone .....	81
Mixing pad .....	79
MKZ EM activator .....	60
MKZ primer .....	60

## N

neo.lign teeth .....	52-55
Living mould chart A and P .....	57
shape kit - reduced quantity of shapes .....	57
novo.lign veneers .....	38-41
Living mould chart novo.lign - novo.lign A & novo.lign P .....	45
shape kit - reduced quantity of shapes .....	45

**O**

Opaque	
crea.lign opaque .....	63
crea.lign opaque Z, zirconium liner .....	63
Opaque combo.lign catalyst.....	62
Opaque combo.lign .....	62

**P**

Polishing finishing set for composite veneers visio.lign toolkit.....	25, 78
polishing paste	
Abraso-Star high lustre polishing .....	78
Acrypol polishing paste for veneering materials .....	78
Primer	
bond.lign stainless steel tripod.....	81
K primer.....	61
MKZ EM activator.....	60
MKZ primer.....	60

**R**

Red-white aesthetics set .....	24
refill pack double-mix cartridge	
DTK adhesive opaque .....	66
DTK adhesive transparent.....	66
Rotary piston crea.lign paste .....	80

**S**

Set adhesive bond	
Chipping repair kit.....	67
combo.lign complete kit .....	66
crea.lign opaque set.....	66
DTK adhesive abutment assortment.....	66
DTK adhesive assortment transparent.....	66
DTK adhesive opaque abutment set.....	66
DTK adhesive opaque assortment .....	66
Full range bonding kit.....	67
refill pack double-mix cartridge DTK adhesive opaque.....	66
refill pack double-mix cartridge DTK adhesive transparent .....	66
Set light-curing device bre.Lux PowerUnit 2 .....	75
Set teeth	
Living mould chart neo.lign - neo.lign A and neo.lign P .....	57
neo.lign shape kit - reduced quantity of shapes .....	57
Set veneer composites	
crea.lign assortment with syringe holder.....	22
crea.lign freestyle intro kit.....	23
crea.lign paste complete kit .....	23
crea.lign paste starter kit, 10 x crea.lign paste.....	22
crea.lign paste starter kit, 9 pastes can be selected.....	22
crea.lign start set N .....	24
crea.lign test set colours A2, A3, A3,5 .....	25
crea.lign veneering system paste + gel, complete set .....	21
Red-white aesthetics set .....	24
visio.lign toolkit finishing set for composite veneers.....	25, 78
visio.paint stains assortment .....	18, 24
Set veneers	
Living mould chart novo.lign - novo.lign A & novo.lign P .....	45
novo.lign shape kit - reduced quantity of shapes .....	45
visio.lign starter kit 4 .....	44
visio.lign veneering system novo.lign A + P starter kit 2 .....	43
visio.lign veneering system starter kit 3 .....	43



<b>Silicone</b>	
Mixing cannulas for visio.sil and visio.sil ILT silicones .....	81
Mixing cannulas for visio.sil fix silicone.....	81
visio.sil fix .....	8
visio.sil ILT .....	9
visio.sil silicone dispenser Multisil dosing device .....	81
visio.sil.....	9
stains visio.paint assortment.....	18, 24
syringe holder for crea.lign paste .....	80
<b>T</b>	
Teeth neo.lign .....	52-55
Thermo-Pen .....	80
Tool stand and glass jar, without tools.....	78
<b>V</b>	
<b>veneer composites</b>	
Application cannulas + pin cover .....	79
crea.lign composite shaper kit .....	80
crea.lign gel dentine .....	16
crea.lign gel GUM .....	16
crea.lign gel Incisal.....	16
crea.lign gel Modifier.....	16
crea.lign gel Transpa Clear.....	16
crea.lign gel .....	12, 16
crea.lign paste dentine .....	17
crea.lign paste GUM.....	17
crea.lign paste Incisal.....	17
crea.lign paste Transpa opal.....	17
crea.lign paste.....	12, 17
crea.lign stainless steel syringe holder .....	73
crea.lign Stains.....	16
crea.lign syringe holder for crea.lign 5 g and stains .....	79
crea.lign syringe holder stainless steel .....	79
Mixing pad.....	79
Rotary piston crea.lign paste.....	80
syringe holder for crea.lign paste .....	80
<b>Veneers</b>	
novo.lign veneers .....	38-41
Thermo-Pen .....	80
visio.CAM composite blanks breCAM.HIPC.....	70
<b>visio.lign</b>	
Living mould chart novo.lign - novo.lign A & novo.lign P .....	45
starter kit 4 .....	44
Toolkit finishing set for composite veneers.....	25, 78
veneering system novo.lign A + P starter kit 2 .....	43
veneering system starter kit 3 .....	43
visio.link.....	61
visio.paint stains assortment.....	18, 24
visio.sil fix .....	8
visio.sil ILT .....	9
visio.sil .....	9







**visio.lign**

The Aesthetic and Functional System

veneering composite



**crea.lign**

veneers



**novo.lign**

denture teeth



**neo.lign**

composite blanks



**visio.CAM**

stains



**visio.paint**

bonding system



**bond.lign**

**bredent**  
group

009504GB-20241129  
Mistake and subject to change reserved

