

breCAM.splint

SPLINT MATERIALS

Occlusal splints for any requirement

For laboratories that rely on precision and clarity

breCAM.splint SOLID

breCAM.splint

breCAM.splint FLEX



Three splint materials – what you gain at a glance:

- **Optimal patient care**
Perfect fit and comfort – for satisfied patients and lasting results.
- **High-quality materials**
Durable, stable and reliable – developed for daily use in the laboratory.
- **Safe & efficient**
Smooth processing with common CAD/CAM systems – for maximum productivity.
- **Biocompatible & well tolerated**
Free of bisphenol A and other questionable ingredients – ideal for sensitive patients.

bredent

Properties & benefits

breCAM.splint SOLID

- + Highly cost-effective
- + Solid

breCAM.splint

- + Optimal cost-benefit ratio
- + Longevity

breCAM.splint FLEX

- + Thermoactive
- + High wearing comfort

Product features			
Wearing time	max. 12 months	max. 2 years	max. 12 months
Biocompatible	Yes	Yes	Yes
Bisphenol A-free	Yes	Yes	Yes

Mechanical properties			
Elasticity modulus	> 1900 MPa	> 2200 MPa	> 400 MPa
Flexural strength	> 60 MPa	> 100 MPa	~ 20 MPa
Water solubility	≤ 0.15 µg/mm ³	≤ 4.5 µg/mm ³	≤ 4.5 µg/mm ³
Water absorption	≤ 25 µg/mm ³	≤ 25 µg/mm ³	≤ 32 µg/mm ³

The right solution for any indication
– precisely tailored to the requirements of dental laboratories.

breCAM.splint SOLID				
REF	Diameter and height	Colour	Material	CE
SPLTSD15	Ø 98.4 x 15 mm	clear	PMMA	CE
SPLTSD20	Ø 98.4 x 20 mm	clear	PMMA	CE

breCAM.splint FLEX				
REF	Diameter and height	Colour	Material	CE
SPLTFX20	Ø 98.4 x 20 mm	clear	Polyacrylate	CE

breCAM.splint				
REF	Diameter and height	Colour	Material	CE
54002315	Ø 98.4 x 15 mm	clear	PMMA	CE 0483
54002320	Ø 98.4 x 20 mm	clear	PMMA	CE 0483
54002325	Ø 98.4 x 25 mm	clear	PMMA	CE 0483
54002516	Ø 95 x 16 mm ZZ	clear	PMMA	CE 0483
54002520	Ø 95 x 20 mm ZZ	clear	PMMA	CE 0483

0091239GB-20250930 Mistake and subject to change reserved



bredent

bredent
group