



# visio.lign

The aesthetic, functional system

## Luting SYSTEM

### The Luting SYSTEM - your smart solution for maximum efficiency in the laboratory



### design natural beauty

#### **digital.connect | digital.bond**

The Luting SYSTEM, consisting of digital.connect and digital.bond, has been developed as a two-component fixation system specially designed for the long-term and secure bonding of prefabricated or CAD/CAM-manufactured prosthetic teeth and tooth segments with the denture base.

powered by  
**visio.lign**

Strong connections - your advantages at a glance:

- **Reliable adhesive bond**  
Enables a stable, long-lasting bond between PMMA dental materials (e.g. breCAM.multiCOM+, neo.lign) and the base (e.g. breCAM.base)
- **Increased efficiency due to elimination of traditional polymerisation**  
No need for hot or cold polymerisation  
- this saves time, energy and material costs
- **Versatile compatibility with modern materials**  
digital.connect & digital.bond are compatible with printed and milled common PMMA and composite materials such as breCAM.multiCOM+ and breCAM.base
- **Clinically tested**  
As a certified class IIa medical device, the system meets high biocompatibility and product safety requirements

# Luting SYSTEM

## Find out more now. Bond now.

### digital.connect

In modern digital dental technology, every detail counts - especially when it comes to bonding. **digital.connect** is the light-curing special primer for a **durable stable bond** between **PMMA bases and prosthetic teeth**.

### What you gain at a glance:

- Reactive boundary layer** - for secure bonding with digital prostheses
- Light curing** - for quick, easy and controlled curing
- Efficient** - 10 ml glass bottle for multiple applications



### digital.bond

With **digital.bond**, you can **position prosthetic teeth in the cavity precisely and without any gaps** - quickly, cleanly and material-safely.

### What you gain at a glance:

- Secure bonding** - between PMMA bases (e.g. breCAM.base) and prosthetic teeth
- Universal colour matching** - pink coloured, perfectly matched to uni.lign PC 20, but also discreetly compatible with other prosthesis colours
- Clean processing** - thanks to practical 8 g double-mix cartridge
- Efficient** - one cartridge is sufficient for approx. 5 prostheses
- Time saving** - dual-curing: self-curing & light-curing



### Luting SYSTEM

digital.connect and digital.bond are also available as a set



Your connection starts here - Ordering information at a glance

REF	Description
VLDDLSET	<b>Luting SYSTEM</b> (contents: 1 x digital.connect, 10 ml glass bottle 1 x digital.bond, 8 g double cartridge 1 x plunger, black 5 x mixing cannula, size 2)

The set is supplied in a sliding box.

REF	Description
VLDDCONN	<b>digital.connect</b> , 10 ml
VLDBBOND	<b>digital.bond</b> , 8 g (contents: 1 x plunger black, 5 x mixing cannula, size 2)

Both products can be ordered individually packaged in a folding box.



0091246GB-20250808  
Mistake and subject to change reserved

bre<sup>dent</sup>

bre<sup>dent</sup>