



visio.lign

The aesthetic, functional system

Luting SYSTEM

The Luting SYSTEM - your smart solution for maximum efficiency in the laboratory



design natural beauty

digital.connect | digital.bond

The Luting SYSTEM, consisting of digital.connect and digital.bond, has been developed as a two-component fixation system specially designed for the long-term and secure bonding of prefabricated or CAD/CAM-manufactured prosthetic teeth and tooth segments with the denture base.

powered by
visio.lign

Strong connections - your advantages at a glance:

- **Reliable adhesive bond**
Enables a stable, long-lasting bond between PMMA dental materials (e.g. breCAM.multiCOM⁺, neo.lign) and the base (e.g. breCAM.base)
- **Increased efficiency due to elimination of traditional polymerisation**
No need for hot or cold polymerisation
- this saves time, energy and material costs
- **Versatile compatibility with modern materials**
digital.connect & digital.bond are compatible with printed and milled common PMMA and composite materials such as breCAM.multiCOM⁺ and breCAM.base
- **Clinically tested**
As a certified class IIa medical device, the system meets high biocompatibility and product safety requirements

Luting SYSTEM

Find out more now. Bond now.

digital.connect

In modern digital dental technology, every detail counts - especially when it comes to bonding. **digital.connect** is the light-curing special primer for a **durable stable bond** between **PMMA bases and prosthetic teeth**.

What you gain at a glance:

- Reactive boundary layer** - for secure bonding with digital prostheses
- Light curing** - for quick, easy and controlled curing
- Efficient** - 10 ml glass bottle for multiple applications



digital.bond

With **digital.bond**, you can **position prosthetic teeth in the cavity precisely and without any gaps** - quickly, cleanly and material-safely.

What you gain at a glance:

- Secure bonding** - between PMMA bases (e.g. breCAM.base) and prosthetic teeth
- Universal colour matching** - pink coloured, perfectly matched to uni.lign PC 20, but also discreetly compatible with other prosthesis colours
- Clean processing** - thanks to practical 8 g double-mix cartridge
- Efficient** - one cartridge is sufficient for approx. 5 prostheses
- Time saving** - dual-curing: self-curing & light-curing



Luting SYSTEM

digital.connect and digital.bond are also available as a set



Your connection starts here - Ordering information at a glance

REF	Description
VLDDLSET	Luting SYSTEM (contents: 1 x digital.connect, 10 ml glass bottle 1 x digital.bond, 8 g double cartridge 1 x plunger, black 5 x mixing cannula, size 2)

The set is supplied in a sliding box.

REF	Description
VLDDCONN	digital.connect , 10 ml
VLDDBOND	digital.bond , 8 g (contents: 1 x plunger black, 5 x mixing cannula, size 2)

Both products can be ordered individually packaged in a folding box.



0091246GB-20250808
Mistake and subject to change reserved

bre^{dent}

bre^{dent}