


JUST PRINT IT.
01 ■ MASTER MODEL



brePRINT

bredent
group

PRECISE SURFACES

01 ■ MASTER MODEL

Are you looking for precision, scratch resistance and high solidity together with high flexural strength? **01 Master Model** guarantees high processing safety in your laboratory as well as excellent results for your customers. Master Resin of the brePRINT product range provides these properties without compromise and offers the highest security during its application. Additionally, the material presents very low shrinkage and an insulating surface against PMMA.



Scan the QR code for the correct print parameters.

■ Properties



Precise ++



Flexural strength



Self-insulating surfaces



Heat resistant ++



Fracture resistant +

■ Scope of application

- Die models
- Models for metal frameworks
- Implant models
- Opposing models
- Hollow models
- Splint models

■ Colours

SAND

LIGHT GREY

■ REF

■ 52000340

■ 52000341

■ Contents

1000 g

0091165GB-20260504
Mistake and subject to change reserved



bre^{dent}

bre^{dent} group