


JUST PRINT IT.
03 ■ DENTAL CAST



brePRINT



3
brePRINT


3D Printing
Resin | DLP

03 ■
DENTAL
CAST

bredent
group

PERFECT FIT

03 ■ DENTAL CAST



Say good-bye to residues when burning-out. 03 Dental Cast of our brePRINT product line is the reliable laboratory partner for residue-free burning-out and presents perfect fitting accuracy. Thanks to the easy handling of Cast Resin of the brePRINT product line your daily amount of laboratory work will be significantly reduced. In addition the material has low shrinkage.



Scan the QR code for the correct print parameters.

■ Properties



Precise +



Smooth surfaces



Residue-free burning out



Self-insulating surfaces



Very low shrinkage

■ Scope of application

- Partial frameworks
- Press ceramic
- Crowns
- Bridges
- Inlays
- Onlays

■ Colours

PURPLE

■ REF

■ 52000360

■ Contents

1000 g

0091167GB-20260504
Mistake and subject to change reserved



bre^{dent}

bre^{dent} group