breCAM.splint FLEX



Perfect fit thanks to
individual adjustment.

Flexible for more comfort.

Flexible. Comfortable. Individual.

It's this simple:

1 Heat.

Before inserting the **breCAM. splint FLEX splint**, immerse it in warm water (**approx. 40°C) for about 1 minute** until it becomes flexible.

2 Adjust.

The heat makes the material flexible so that it adapts optimally to your dental situation.

3 Wear.

Once cooled, the customised shape is retained and ensures that it is comfortable to wear.

bredent group

0091198GB_20260925